


**HEALTHMARK ECONOMY TRAYS****Product Code(s)**

- HM30732-H
- HM30842-H
- HM31042-H
- HM31462-H
- HM31772-H



A range of hinged plastic trays constructed with an economy budget in mind.

The perforations are designed in such a way that the inserts from the PST line will fit neatly in each tray. These trays are ideal when it is not necessary to view the contents through the lid.

All trays are tan coloured except for the HM31772-H which is blue.

**Material**

All Healthmark Economy Trays are manufactured from Polypropylene.

**Sizes**

HM30732 – H	-184 x 83 x 29mm
HM30842 – H	-222 x 114 x 44mm
HM31042 – H	-260 x 114 x 44mm
HM31462 – H	-356 x 159 x 64mm
HM31772 – H	-445 x 197 x 60mm

**Sterilising**

Healthmark Economy Trays can be steam sterilised at temperatures up to 134°C. They can also be sterilised by ethylene oxide and plasma. They can also be wash-decontaminated.