

PACK AND TRAY STRAPS

Product Code(s)

- STR002
- STR004



An easy solution to the problems of packs blowing during autoclaving. These straps are extremely strong with double welded joints. They also offer a safe means of lifting heavy trays without risk of damage to wraps. The straps are held in place with hook and loop fastenings.

Material

Pack and Tray Straps are manufactured from nylon.

Colour and Sizes

STR002	Yellow	30 x 650mm
STR004	Blue	30 x 1400mm

Sterilisation

Pack and Tray Straps are suitable for steam autoclaving up to 137°C, ethylene oxide and plasma sterilisation processes up to 60°C.

Each strap has a 10 cm "hook area" to bind to the Velcro strap.