

TRAY TAGS

Product Code(s)

- TAG011 – TAG027 inclusive (Pack of 100 Tags)
- TAG011-TAG027/P20 inclusive (Pack of 20 Tags)
- TAG55WH
- TAGMULTIPK1
- TAGMET7530

Coloured plastic labelling tags (except TAGMET7530 which is made of stainless-steel) for tracking and tracing purposes.

Split rings and plastic security Ties are available for attaching Tray Tags to items such as stainless-steel baskets.

Plastic tray tags are available in 17 different colours.

Material

Plastic Tray Tags by an injection moulding process using polymer Polyamide 6.

TAGMET7530 is made of AISI304 Non-Magnetic Stainless Steel.



All tags are latex-free.

Sizes

TAG011 – TAG027	90mm x 38mm
TAG55WH	90mm x 55mm
TAGMET7530	90mm x 30mm

Specification**Plastic:**

Glass Content	30 %
Density	1.34-1.36 g/cm ³
Water Absorption	6.0 %
Tensile Strength	165 N/mm ²
Elongation at Break	2.8 %
Melt Temp.	222 °C

Stainless Steel:

Density	8 kg/dm ³
Heat Capacity	500 J/kg°C
Thermal Conductivity	15 W/m°C
Electrical Resistivity	0.75μΩm

Sterilisation compatibility

Plastic Tray Tags are suitable for use in steam (up to 137°C) sterilisation processes. The Tags will withstand multiple sterilisation cycles- in particular they will withstand at least 100 steam sterilisation cycles.

Plastic Tray Tags can be used in plasma sterilisation processes however the life of the Tray Tag will be significantly reduced compared to use within steam or ethylene oxide sterilisation processes.

TAGMET7530 is compatible to steam, Ethylene Oxide and Plasma sterilisation methods.

TRAY TAG COLOURS

CODE	PRODUCT
TAG011	White
TAG012	Tan
TAG013	Yellow
TAG014	Orange
TAG015	Red
TAG016	Blue
TAG017	Green
TAG018	Grey
TAG019	Purple
TAG020	Black
TAG021	Lime
TAG022	Pink
TAG023	Turquoise
TAG024	Mustard
TAG025	Terracotta
TAG026	Royal Blue
TAG027	Sky Blue
TAG55WH	White