

TRANSPORTATION TRAYS**Product Code**

- TTRAY001
- TTRAY002
- TTRAY003

A range of trays designed to protect wrapped packs from external damage and to reduce damage from stacking.
May used alone or with a strap.
Can be used either singly, or as a pair to form a lidded box.

**Material**

TTRAY001 and TTRAY003 are manufactured from cardboard (cellulose fibre).
TTRAY002 is manufactured from Polypropylene.

Dimensions

- | | |
|------------|------------------|
| ➤ TTRAY001 | 500 x 270 x 70mm |
| ➤ TTRAY002 | 500 x 270 x 70mm |
| ➤ TTRAY003 | 500 x 540 x 70mm |

Cleaning and Sterilisation

TTRAY002 may be wash decontaminated in any wash decontaminating equipment.

TTRAY001 and TTRAY003 are NOT suitable for wash decontamination.

TTRAYs are NOT suitable for sterilising.