

CLINIPAK TRANSPORTATION COVER BAGS

Product Code(s)

- TCOVER001
- TCOVER002
- TCOVER003
- TCOVER007
- TCOVER008
- TCOVER-RED001
- TCOVER-RED002
- TCOVER-RED003
- TCOVER-YELL003



Clinipak Transportation Cover Wrap available in bag format in a variety of colours into which barrier-wrapped trays are inserted. The Bag can be closed with the lip provided and closed with autoclave tape. The Clinipak Transportation Cover Wrap is available in the colours Purple, Red and Yellow. The material is inherently water repellent and fully steam permeable, with a permeability level 400% higher than a typical nonwoven barrier fabric. **Clinipak Transportation Cover Bags do not contain a dedicated microbiological barrier layer and MUST always be used in combination with barrier wrap(s) conforming to EN 868 part 1 and ISO 11607.**

Material

Clinipak Transportation Cover Bags are made from 100% Polypropylene.

Sizes

- | | |
|---|--|
| ➤ TCOVER001 and TCOVER-RED001 | 400X700mm + 90mm lip - 200 bags per pack |
| ➤ TCOVER002 and TCOVER-RED002 | 500x700mm + 90mm lip - 200 bags per pack |
| ➤ TCOVER003, TCOVER-RED003 and TCOVER-YELL003 | 600x700mm + 90mm lip - 100 bags per pack |
| ➤ TCOVER007 | 400x550mm + 90mm lip - 100 bags per pack |
| ➤ TCOVER008 | 550x800mm + 90mm lip - 100 bags per pack |

Sterilising

Clinipak Transportation Cover Bags are compatible with steam sterilisation processes up to 137°C,
Clinipak Transportation Cover Bags are also compatible with both Ethylene Oxide and Plasma sterilisation systems.

Standard Properties

<i>Basis Weight</i>	<i>60 gsm</i>
<i>MD Tear</i>	<i>7 N</i>
<i>Tensile Strength MD (dry)</i>	<i>1.65 kN/m</i>
<i>Tensile Strength CD (dry)</i>	<i>1.13 kN/m</i>
<i>Elongation MD</i>	<i>115%</i>
<i>Elongation CD</i>	<i>115%</i>
<i>Microbiological Barrier (dry)</i>	<i>No Barrier</i>
<i>Microbiological barrier (wet)</i>	<i>No Barrier</i>