

CLINIPAK TRANSPORTATION COVER DRY BAGS

Product Code(s)

- TCOVERDRY001
- TCOVERDRY002
- TCOVERDRY003

Clinipak Transportation Cover Dry Wrap available in bag format into which barrier-wrapped trays are inserted. The Bag can be closed with the lip provided and closed with autoclave tape. The absorbent surface of the wrap promotes dry packs. The Clinipak Transportation Cover Dry Wrap is only available in Purple colour.



The material is fully steam permeable with a permeability level 300% higher than a typical nonwoven barrier fabric. **Clinipak Transportation Cover Dry Bags do not contain a dedicated microbiological barrier layer and MUST always be used in combination with barrier wrap(s) conforming to EN 868 part 1.**

Material

Purple layer made from 100% polypropylene.

White layer made from Polypropylene and viscose.

Sizes

- | | | |
|----------------|----------------------|-------------------|
| ➤ TCOVERDRY001 | 400X700mm + 90mm lip | 200 bags per pack |
| ➤ TCOVERDRY002 | 500x700mm + 90mm lip | 200 bags per pack |
| ➤ TCOVERDRY003 | 600x700mm + 90mm lip | 100 bags per pack |

Sterilising

Clinipak Transportation Cover Dry Bags are steam sterilisable.

Clinipak Transportation Cover Dry Bags are also compatible with Ethylene Oxide sterilisation.

Standard Properties

Basis Weight	75 gsm
Tensile Strength MD (dry)	1.8kN/m
Tensile Strength MD (wet)	0.8kN/m
Elongation MD	60%
Elongation CD	80%
Surface Absorbency	85g/m ² (white surface)
Permeability	500 l/min/100cm ²
MD Tear	10N
Microbiological Barrier (dry)	No Barrier
Microbiological barrier (wet)	No Barrier