

**CLINIPAK PURPLE MATS****Product Code(s)**

➤ MAT001 – MAT999

Extremely strong Clinipak Purple Mats designed to sit underneath trays and baskets for added protection against sharp-edged trays.

**Material**

Clinipak Purple Mats are made from 100% Polypropylene.

**Sizes**

A wide variety of sizes are available to fit the most common sizes of trays.

**Sterilising**

Clinipak Purple Mats are suitable for use in steam sterilisation up to 137°C, ethylene oxide, and Plasma (Hydrogen peroxide) sterilising systems. (e.g. Steris V Pro)

**Compliance**

This material does not comply with EN 868 Parts 1 & 2 or ISO 11607-1.

**Standard Properties**

Basis Weight	60 gsm
MD Tear	7 N
Tensile Strength MD (dry)	1.65 kN/m
Tensile Strength CD (dry)	1.13 kN/m
Elongation MD	115 %
Elongation CD	115 %
Microbiological Barrier (dry)	No Barrier
Microbiological barrier (wet)	No Barrier

**Storage/shelf life**

Purple Mats have a shelf life of 5 years prior to use when stored at 20-80%RH and 10-35°C.