

CLINIPAK PURPLE DRY MATS

Product Code(s)


➤ DRYMAT001 –DRYMAT999

Strong spun polypropylene sheet with an absorbent layer for use as a strong mat under trays - to protect against sharp-edged trays and provide rapid moisture dispersion during steam sterilisation.

Most commonly used with a white crepe tray liner inside the tray.

Material

Clinipak Purple Dry Mats are made from polypropylene and viscose.



Sizes

A wide variety of sizes are available to fit the most common sizes of trays.

Sterilising

Clinipak Purple Dry Mats are designed primarily for use in steam sterilisation up to 137°C.

Clinipak Purple Dry Mats are not suitable for use in Plasma sterilising systems or ethylene oxide sterilisation systems.

Standard Properties

Basis Weight	83 gsm max.
Tensile Strength MD	180 N
Tensile Strength CD	80 N
Elongation MD	70 %
Elongation CD	120%
MD Tear	3.7 N min.
Microbiological Barrier (dry)	No Barrier
Microbiological barrier (wet)	No Barrier

Storage/Shelf Life

Clinipak Purple Dry Mats have a shelf life of 5 years prior to use when stored at 20-80%RH and 10-35°C.