

POUCH AND BAG RACK**Product Code(s)**

- DISP002

Specially designed racks for neatly holding pouches and bags on a working surface.

**Contents**

A DISP002 pouch and Bag Rack includes the following items:

- Plastic Dispenser container 3 large racks and 1 small rack
- Black PEN001
- Red PEN002
- Blue PEN003
- Pen(s) holder

Materials

Pouch and Bag Rack – Perspex {Poly(methyl methacrylate)}

Dimensions

DISP002 large rack size = 635 x 315 x 335mm
 Small Rack size =630 x 120 x 85mm