

### Product Codes

- CBL000 – CBL999
- CBL000X000 – CBL999X999

Clinipak Cut Corner Tray Liners are highly absorbent crepe tray liners produced from pure cellulose wood pulp, manufactured with a low chloride and sulphate content and free from optical dyes. Clinipak Corner Cut Tray Liners have cuts at the corners to allow the edges to be folded up. The liners can be safely used inside trays, next to delicate instruments for steam sterilization or may be used underneath trays. In both cases, condensate dispersion facilitates the production of dry packs. Clinipak Cut Corner Tray Liners are designed for use in steam sterilising systems, although they are also compatible with ethylene oxide sterilising systems. They are not suitable for use in plasma sterilising systems.

### Standard Properties

Test method			Measure	min	max
Substance	DIN 53104	-	g/m <sup>2</sup>	78	82
Thickness	DIN 53105 T1	-	µm	-	-
Tensile strength	DIN 53112 T1	md	N/15mm	35,3	-
Tensile strength	DIN 53112 T1	cd	N/15mm	22,4	-
Wet strength	DIN 53112 T2	md	N/15mm	10,0	-
Wet strength	DIN 53112 T2	cd	N/15mm	6,0	-
Breaking length	DIN 53112 T1	md	km	3,0	-
Breaking length	DIN 53112 T1	cd	km	1,9	-
Stretch	DIN 53112 T1	md	%	8	-
pH	online	-	-	6	8
Ash content	DIN 54370	-	%	0	1
Absolute moisture content	online	-	%	3,5	6,5
Capillary rise	DIN 53106	md	mm	40	-
Capillary rise	DIN 53106	cd	mm	45	-
Brightness (white Paper)	DIN 53145 - R457	-	%	82	90

ISO 9001; ISO 14001; ISO 50001; REACH

MD=Machine Direction,CD=Cross Direction,T=Top,B=Bottom

All above are typical values at conditioning atmosphere 50% and 23°C according to ISO-standards.

### Sterilising Compatibility

Clinipak Cut Corner Tray Liners are suitable for use in steam sterilising systems, up to 137°C/4min, and ethylene oxide sterilising systems.

They are not suitable for use in plasma sterilising systems.

**CLINIPAK CUT CORNER TRAY LINER DESIGNS**

Product Code	Overall Size (mm)	Size Cut "A"	Size Cut "B"
CBL003	13.5 x 23"	60mm	60mm
CBL300X300	300 x 300mm	50mm	50mm
CBL350X350	350 x 350mm	60mm	60mm
CBL400X400	400 x 400mm	60mm	60mm
CBL400X600	400 x 600mm	60mm	60mm
CBL500X700	500 x 700mm	60mm	60mm
CBL580X350	580 x 350mm	60mm	60mm
CBL580X350H	580 x 350mm	50mm	50mm

