

AUTOCLAVE TAPE DISPENSERS**Product Code(s)**

- DISP004 3-Reel Dispenser
- DISP005 2-Reel Dispenser

Tape dispensers for easy application of autoclave tape.
An integral cutter-guard is included for added safety.

**Materials**

Autoclave Tape Dispenser – metal and plastic

Dimensions and weights

DISP004 - 85(W) x 230(L) x 105(H) mm weight = 1.40kg

DISP005 - 80(W) x 230(L) x 105(H) mm weight = 1.30kg