

SINGLE-USE AUTOCLAVE TAGS

Product Code(s)

- TAG-AUTO1
- TAG-AUTO1GRN, BLU, RED, ORA, YEL

Single-Use Autoclave Tags are designed for the identification of trays that are being steam sterilised.



Material

Single-Use Autoclave Tags are manufactured from polyester.

Properties

Size	50 x 228.6 mm
Thickness	175 micron
Weight	245 gsm
Colours	White, green, blue, red, orange and yellow

Printability

Variable information printing and barcodes can be put on both sides using thermal transfer printers. A Black PEN001 may also be used to hand-write information on the tag.

Pre-printed Tags

As an alternative to on-demand printing of blank tags, Single-Use Autoclave Tags can be supplied pre-printed.

Sterilisation Compatibility

Single-Use Autoclave Tags are compatible with steam sterilisation systems up to 137degC/3.5min.
Single-Use Autoclave Tags are not compatible with other sterilisation systems.