

SINGLE-USE CONTAINER CARRIER TAGS

Product Code(s)

- TAG-CONTAIN1
- TAG-CONTAIN2
- TAG-CONTAIN1-GRE
- TAG-CONTAIN1-RED
- TAG-CONTAIN-YEL

Single-Use Container Carrier Tags are designed for the identification of containers that are to be steam sterilised. Two sizes are available (see below for details), both sizes are available in white, in addition the size TAG-CONTAIN1 is also available in green, red and yellow colours.



Material

Single-Use Container Carrier Tags are manufactured from polyester.

Properties

Size	TAG-CONTAIN1xxx	75 x 228.6 mm,
	TAG-CONTAIN2	96 x 280 mm
Thickness		175 micron
Weight		245 gsm
Colours		white, green, red and yellow

Printability

Variable information printing and barcodes can be put on both sides using thermal transfer printers. A black PEN001 may also be used to hand-write information on the Tag.

Pre-printed Tags

As an alternative to on-demand printing of blank tags, Single-Use Container Carrier Tags can be supplied pre-printed.

Sterilisation Compatibility

Single-Use Container Carrier Tags are compatible with steam sterilisation systems up to 137degC/3.5min. Single-Use Container Carrier Tags are not compatible with other sterilisation systems.