

SINGLE-USE SECURITY TAGS**Product Code(s)**

- TAG-SECURE1
- TAG-SECURE1GRN, BLU, RED, ORA, YEL

Single-Use Security Tags are designed for the identification of containers that are to be steam sterilised. Their design enables them to be securely fixed to a container via the standard container hinged lock.

**Material**

Single-Use Security Tags are manufactured from polyester.

Properties

Size	152 x 101.6 mm
Thickness	175 micron
Weight	245 gsm
Colours	white, green, blue, red, orange and yellow

Printability

Variable information printing and barcodes can be put on both sides using thermal transfer printers. A black PEN001 may also be used to hand-write information on the tag.

Pre-Printed Tags

As an alternative to on-demand printing of blank tags, Single-Use Security Tags can be supplied pre-printed.

Sterilisation Compatibility

Single-Use Security Tags are compatible with steam sterilisation systems up to 137degC/3.5min. Single-Use Security Tags are not compatible with other sterilisation systems.