

## SINGLE-USE AUTOCLAVE TRAY TAGS

### Product Code(s)

- TAG-TRAY1
- TAG-TRAY1GRN, BLU, RED, ORA, YEL

Single-Use Autoclave Tray Tags are designed for the identification of trays that are to be steam sterilised.



### Material

Single-Use Autoclave Tray Tags are manufactured from polyester.

### Properties

Size	70 x 101.6 mm
Thickness	175 micron
Weight	245 gsm
Colours	white, green, blue, red, orange and yellow

### Printability

Variable information printing and barcodes can be put on both sides using thermal transfer printers. A black PEN001 may also be used to hand-write information on the tag.

### Pre-Printed Tags

As an alternative to on-demand printing of blank tags, Single-Use Autoclave Tray Tags can be supplied pre-printed.

### Sterilisation Compatibility

Single-Use Autoclave Tray Tags are compatible with steam sterilisation systems up to 137degC/3.5min. Single-Use Autoclave Tray Tags are not compatible with other sterilisation systems.