


**SINGLE USE SELF-TIE TAGS****Product Code(s)**

- SELFTIE01
- SELFTIE01-GRN, BLU, RED, ORA, YEL & USED

Single-Use Self-Tie Tags are designed for the identification of trays that are to be steam sterilised.

**Material**

Single-Use Self-Tie Tags are manufactured from polyester.

**Properties**

Size	25 x 190.5 mm
Thickness	175 micron
Weight	245 gsm
Colours	white, green, blue, red, orange and yellow

**Printability**

Variable information printing and barcodes can be put on both sides using thermal transfer printers. A black PEN001 may also be used to hand-write information on the Tag.

**Pre-Printed Tags**

As an alternative to on-demand printing of blank tags, Single-Use Self-Tie Tags can be supplied pre-printed.

**Sterilisation Compatibility**

Single-Use Self-Tie Tags are compatible with steam sterilisation systems up to 137degC/3.5min.  
Single-Use Self-Tie Tags are not compatible with other sterilisation systems.