

CBM STAINLESS STEEL BASKETS**Product Code(s)**

➤ CBMxxxx

CBM Stainless Steel Baskets are range of sterilisation baskets designed for the washing and sterilisation of small and very small items. The baskets are generally supplied with a lid which incorporates a sliding catch. The lid is laser-etched with LOT number.

A range of accessories are available, including silicone instrument holders.

**Material**

CBM Stainless Steel Baskets are manufactured from electro-polished Grade 304 Stainless Steel.

Decontamination and sterilisation

CBM Stainless Steel Baskets can be disinfected using standard equipment and conditions such as washing in a deluge washing machine at 90°C using an enzymatic detergent.

CBM Stainless Steel Baskets can be sterilised using steam in a porous load autoclave at 134 - 137°C for 3-3.5 minutes.