

**CHUX TURQUOISE CREPED LOW LINTING WIPES****Product Code(s)**

➤ CHUX-WIPE01

**Product Description/Intended uses**

CHUX Turquoise Creped Low Linting Wipes are used in the wash room to remove moisture from manually-cleaned instruments /medical devices. These are also used in the inspection, assembly and packing room (clean room), to remove excess moisture droplets from instruments following automated cleaning and disinfection. Their high water-absorbing and low-linting properties ensure that instruments are dry and fibre-free.

**Material**

CHUX-WIPE01 wipes are manufactured from 56%Pulp/44%Polyester

**Properties-Typical values**

|                           |                     |
|---------------------------|---------------------|
| Weight                    | 70 g/m <sup>2</sup> |
| Thickness                 | 2.1 mm/4plies       |
| Sink time                 | <2sec               |
| Water Absorption Capacity | 900%                |

**Size**

|            |               |
|------------|---------------|
| Folded     | - 155 x 166mm |
| Opened out | - 310 x 330mm |

**Single use**

CHUX-WIPE01 wipes can be used to remove moisture from many instruments.

Once they have been used, they should not be re-used.

**Sterilising Compatibility**

CHUX Turquoise Creped Low Linting Wipes are single use and not suitable for sterilisation.