


PRE-PRINTED STERILISABLE SELF-ADHESIVE LABELS

Product Code

- | | |
|------------------|--|
| ➤ SAL-FASTTRACK | - States "FAST TRACK" |
| ➤ SAL-HEAVYSET | - States "HEAVY SET" |
| ➤ SAL-LOANSET | - States "LOAN SET" |
| ➤ SAL-FRAGILE | - States "FRAGILE INSTRUMENTS" |
| ➤ SAL-DISINFECT | - States "DISINFECTED" |
| ➤ SAL-DONOTSTACK | - States "DO NOT STACK" |
| ➤ SAL-MISSING | - States "MISSING INSTRUMENT" |
| ➤ SAL-SERVICE | - States "SERVICE REPORT RAISED DO NOT OPEN" |
| ➤ SAL-URGENT | - States "URGENT" |



Pre-Printed Sterilisable Self-Adhesive Labels for wrapped instrument sets which require special identification in an SSD unit. Designed to be applied to the outer sterilisation wrap, the labels are large and bright, clearly highlighting the label to staff.

Materials

Adhesive	Latex-free self-adhesive formulation
Label stock	White coated matt paper
Release backing paper	Silicone-coated glassine

Label Dimensions

Label size = 76 x 76mm.
1000 labels per roll.

Sterilisation Compatibility

Pre-Printed Sterilisable Self-Adhesive Labels are designed to be used only with steam sterilisation processes up to 137°C/3.5min.

Storage Conditions

Pre-Printed Sterilisable Self-Adhesive Labels should be stored in a cool dry place.