

GREEN CONTAINER PADLOCK (with indicator)**Product Code(s)**

- LOCKGRN605MIND

A plastic container lock with a temperature indicator dot which enables the user to determine whether the container has been through a steam sterilisation process.

**Material**

The lock is manufactured from polypropylene.

Size

420 X 330 X 250mm

Standard Properties

- Latex-free, pvc-free, phthalate-free.
- Indicator lead-free.
- Label stock with process indicator sensitive to colour change.
- Indicator dot changes colour to BLACK.
- Autoclave cycles of (typically) 121°C/20 min or 134° - 137°C/3.5min.
- Inbuilt break points to allow easy removal.
- Can be used in steam sterilisation process.
- Storage best at room temperature with no direct sun or heat exposure.
- Shelf life for indicator locks- three years from manufacture date.

Temperature indication

The indicator colour change occurs at 121°C and higher. The green locks are, therefore, able to indicate that the container has been through a steam sterilisation process.

The locks cannot be used to determine that the steriliser was operating at the correct temperature.

Compliance

Container Padlocks comply with EN 867 Part 1 Class A and with ISO 11140-1.