

**CHUX WHITE LOW LINTING WIPES****Product Code(s)**

➤ CHUX-VERA02A

**Product Description/Intended Use**

CHUX White Low Linting Wipes are used in the wash room to remove moisture from manually-cleaned instruments /medical devices. These are also used in the inspection, assembly and packing room (clean room), to remove excess moisture droplets from instruments following automated cleaning and disinfection.

Their high water-absorbing and low-linting properties ensure that instruments are dry and fibre-free.

CHUX White Low Linting Wipes are SINGLE-USE wipes.

**Material**

CHUX White Low Linting Wipes are manufactured from 56% pulp / 44% polyester.

**Properties - Typical values**

Weight	68 g/m <sup>2</sup>
Thickness	1.6 mm/4plies
Sink time	<2sec
Water Absorption Capacity	755%

**Size**

Folded - 150 x 170mm  
Opened out - 300 x 340mm

**Quantity**

300 pcs per pack.

**Sterilising Compatibility**

CHUX White Low Linting Wipes are not suitable for sterilisation.