

CLINIPAK RELIANCE GRN/BLU WRAP (RELIANCE 340/345)**Product Code(s)**

- RELBLUxxxx
 - RELGRNxxxx
 - RELGRN/BLUxxxx
 - REL340xxxx
- (where xxxx= sheet size)



Clinipak Reliance GRN/BLU Wrap is a latex-free, high-performance sterilisation wrap suitable for inner and outer wrapping of trays and packs. It offers a combination of high strength, good barrier and softness in conjunction with alcohol repellency, antistatic and low linting capabilities.

This material is manufactured from cellulosic and synthetic fibres bonded with an acrylic binder. This structure significantly improves its strength while increasing its softness. It is available in Green and Blue and an interleave of the 2 colours.

Sterilisation

Clinipak Reliance GRN/BLU Wrap is designed for Steam sterilisation at 134-137°C and ETO, Formaldehyde and Irradiation sterilisation methods.

It is NOT suitable for Plasma sterilisation.

Compliance

Clinipak Reliance GRN/BLU Wrap conforms to ISO 11607-1 & BS EN 868-2.

Typical Data

<i>Basis Weight</i>	<i>57gsm. +/- 5%</i>
<i>Extract pH</i>	<i>5.7 pH</i>
<i>Chloride content</i>	<i>0.05% max.</i>
<i>Sulphate content</i>	<i>0.25% max.</i>
<i>Dry Tensile MD</i>	<i>1.65 kN/m</i>
<i>Dry Tensile CD</i>	<i>1.25 kN/m</i>
<i>Elongation MD</i>	<i>8 %</i>
<i>Elongation CD</i>	<i>9 %</i>
<i>Alcohol repellent</i>	<i>Yes</i>

Storage/Shelf Life

Clinipak Reliance GRN/BLU Wrap has a shelf life of 5 years prior to use when stored at 35-50%RH and 15-25°C.