





Product Code(s)

- CC001
- CC002
- H-CC001
- H-CC002



Simple to assemble corner protector for use during sterilisation to eliminate the risk of wrap damage caused by baskets with sharp corners. Assemble by folding and sliding under each corner.

Material

- CC001 / H-CC001 White rigid Card
- CC002 / H-CC002 Soft white crepe paper

Properties**CC001 / H-CC001**

White rigid card is wet strengthened and manufactured with low chloride and sulphate contents and free from optical brightening agents.

Grammage	171 g/m ²
Thickness	220µm
Burst	840 kPa min
Cobb	25 g/m ² max
Wet Burst	250 kPa min
Chloride Content	0.05% max
Sulphate Content	0.25% max
pH of Hot Water extract	pH between 5.0 and 8.0
OBA content	Optical Brightener-Free

CC002 / H-CC002

Soft crepe material is a white absorbent material with low chloride and sulphate contents and free from optical brightening agents.

Grammage	140 g/m ²
Thickness	480 µm
Dry Burst	250 kPa
Wet Burst	90 kPa
Chloride Content	0.023% max
Sulphate Content	0.029% max
Absorbency Capacity	200 g/m ²

Sterilisation Compatibility

Corners are suitable for use in steam (up to 137°C) sterilisation processes. They are not suitable for use in plasma sterilisation processes.