

INSTRUMENT MARKING TAPE

Product Code(s)

- | | |
|-------------------|-----------------------|
| ➤ MKT001 – MKT004 | Ring Binders |
| ➤ MKT005 – MKT015 | Plain colours |
| ➤ MKT016 – MKT025 | Zebra-striped colours |
| ➤ MKT031 | Tiger-striped colours |

Latex-free instrument marking tape is designed to help organise and identify instruments as they move between sterile processing and other departments.

Available in a range of plain and zebra/tiger-striped colours.



Material

Face material: Polymeric stabilized calendared vinyl (PVC), 66 µm thick, cadmium-free.

Adhesive: High quality clear solvent based permanent acrylic adhesive.

Liner: Highly stable 135-gsm white kraft paper with blue back printing.

For white vinyl, the liner is coated with a blue silicone coating to improve the contrast with the film colour.

Properties

Tensile strength: ≥ 24 N/15mm

Elongation at break: ≥ 150 %

Peel adhesion: Peel adhesion 180° on brushed stainless steel (N/25mm)
 - 20 min. residence: 16
 - 24 hours residence: 18
 - 1 week residence: 20

Peel adhesion 180° on mirror stainless steel (N/25mm)
 - 20 min. residence: 11
 - 1-week residence: 13

Solvent and chemical resistance: Resistant to most oils and grease, fuel, aliphatic solvents, weak acids, salts and alkalis.

Temperature ranges: Application: + 8°C to + 40°C
 End-use: - 40°C to + 90°C

Shelf Life

2 years when stored at 15 to 25°C and ± 50 % relative humidity.

Document Ref: CDS-3041

Page 1 of 2

Effective Date: 2020-03-13

TMP-CDS-V02: Effective 12 Feb 2018

Version No: 03

Colours

CODE	PRODUCT
MKT001	Ring Binder with Index only
MKT002	Ring Binder with index and 1 sheet of each plain colour
MKT003	Ring Binder with index and 1 sheet of each zebra colour
MKT004	Ring Binder with index and 1 sheet of each plain and zebra colour
MKT005	Plain Purple
MKT006	Plain White
MKT007	Plain Beige/Tan
MKT008	Plain Yellow
MKT009	Plain Orange
MKT010	Plain Red
MKT011	Plain Blue
MKY012	Plain green
MKT013	Plain Silver
MKT014	Plain Brown
MKT015	Plain Black
MKT016	Zebra Beige/Tan
MKT017	Zebra Yellow
MKT018	Zebra Orange
MKT019	Zebra Red
MKT020	Zebra Blue
MKT021	Zebra Green
MKT022	Zebra Silver
MKT023	Zebra Brown
MKT024	Zebra Black
MKT025	Zebra Purple
MKT031	Tiger Orange

Sterilisation

Instrument Marking Tape is compatible with steam sterilisation up to 137°C for up to 4 min.
Instrument Marking Tape is not compatible with any other type of sterilisation process.