

## TRANSPORTATION BAGS

## Product Code(s)

- TBAG001
- TBAG002
- TBAG003
- TBAG004
- TBAG005
- TBAG006

Polyethylene bags for protecting sterile packs against handling whilst providing an effective dust cover. They have a tamper-proof seal and an easy tear-opening facility.



## Material

Transportation Bags, TBAG001,002,003,004 and 005 are constructed from 70 micron polyethylene.

TBAG006 is constructed from 93 micron polyethylene.

The adhesive is a permanent hot melt adhesive.

## Sizes

TBAG001	380 x 760mm	250 Bags per box
TBAG002	585 x 765mm	200 Bags per box
TBAG003	380 x 710mm	250 Bags per box
TBAG004	385 x 490mm	250 Bags per box
TBAG005	664 x 762mm	200 Bags per box
TBAG006	585 x 765mm	150 Bags per box

## Design

- 1) Permanent hot melt adhesive forms a tamper-proof seal
- 2) With half-moon perforations to allow packs to breathe preventing condensation
- 3) Non-slip polyethylene
- 4) Perforated end to allow easy tear-opening
- 5) Limited Printing to allow content information to be visible

## Sterilising

Transportation bags are not steam sterilisable and should be used only after the pack has been separately steam sterilised.