

Product Codes

- CBL000 – CBL999
- CBL000X000 – CBL999X999

Clinipak Cut Corner Tray Liners are highly absorbent crepe tray liners produced from pure cellulose wood pulp, manufactured with a low chloride and sulphate content and free from optical dyes. Clinipak Corner Cut Tray Liners have cuts at the corners to allow the edges to be folded up. The liners can be safely used inside trays, next to delicate instruments for steam sterilization or may be used underneath trays. In both cases, condensate dispersion facilitates the production of dry packs.

Standard Properties

| Test method | | | Measure | min | max |
|---------------------------------------|------------------|----|------------------|------|-----|
| Substance | DIN 53104 | - | g/m ² | 78 | 82 |
| Thickness | DIN 53105 T1 | - | µm | - | - |
| Tensile strength | DIN 53112 T1 | md | N/15mm | 35,3 | - |
| Tensile strength | DIN 53112 T1 | cd | N/15mm | 22,4 | - |
| Wet strength | DIN 53112 T2 | md | N/15mm | 10,0 | - |
| Wet strength | DIN 53112 T2 | cd | N/15mm | 6,0 | - |
| Breaking length | DIN 53112 T1 | md | km | 3,0 | - |
| Breaking length | DIN 53112 T1 | cd | km | 1,9 | - |
| Stretch | DIN 53112 T1 | md | % | 8 | - |
| pH | online | - | - | 6 | 8 |
| Ash content | DIN 54370 | - | % | 0 | 1 |
| Absolute moisture content | online | - | % | 3,5 | 6,5 |
| Capillary rise | DIN 53106 | md | mm | 40 | - |
| Capillary rise | DIN 53106 | cd | mm | 45 | - |
| Brightness (white Paper) | DIN 53145 - R457 | - | % | 82 | 90 |
| ISO 9001; ISO 14001; ISO 50001; REACH | | | | | |

MD=Machine Direction,CD=Cross Direction,T=Top,B=Bottom

All above are typical values at conditioning atmosphere 50% and 23°C according to ISO-standards.

Sterilising Compatibility

Clinipak Cut Corner Tray Liners are suitable for use in steam sterilising systems, up to 137°C/4min, and ethylene oxide sterilising systems.

CLINIPAK CUT CORNER TRAY LINER DESIGNS

| Product Code | Overall Size (mm) | Size Cut "A" | Size Cut "B" |
|--------------|-------------------|--------------|--------------|
| CBL003 | 13.5 x 23" | 60mm | 60mm |
| CBL300X300 | 300 x 300mm | 50mm | 50mm |
| CBL350X350 | 350 x 350mm | 60mm | 60mm |
| CBL400X400 | 400 x 400mm | 60mm | 60mm |
| CBL400X600 | 400 x 600mm | 60mm | 60mm |
| CBL500X700 | 500 x 700mm | 60mm | 60mm |
| CBL580X350 | 580 x 350mm | 60mm | 60mm |
| CBL580X350H | 580 x 350mm | 50mm | 50mm |

