

CLINIPAK RELIANCE® DEXTEX 200 BLUE/GREEN**Product Code(s)**

- RELIANCE® DEXTEX 200 BLUE (RELBLUxxx)
- RELIANCE® DEXTEX 200 GREEN (RELGRExxxx)
(Where xxxx= sheet size)



Clinipak RELIANCE® DEXTEX 200 BLUE/GREEN is a latex-free, high-performance sterilisation wrap suitable for inner and outer wrapping of trays and packs. It offers a combination of high strength and softness with excellent bacterial barriers in conjunction with alcohol repellency, anti-static and low lint capabilities.

The cellulosic fibers create excellent sterility maintenance where synthetic binders and PET fibers bring more strength and softness to the product.

RELIANCE® DEXTEX 200 is available in two colors - blue and green.

Materials

Manufactured from 70% cellulose, 20% synthetic binder and 10% PET fibers

Sterilisation Compatibility

Clinipak RELIANCE® DEXTEX 200 BLUE/GREEN is compatible with the following sterilisation methods:

- Steam (up to 135°C)
- Ethylene Oxide (EtO)
- Formaldehyde
- Gamma Irradiation (Co 60 - 50 kGy)

It is NOT suitable for plasma sterilisation.

Standard Properties

Basis Weight	57gsm. +/- 5%
Extract pH	5.7 pH
Chloride content	0.05% max.
Sulphate content	0.25% max.
Dry Tensile MD	1.65 kN/m
Dry Tensile CD	1.25 kN/m
Elongation MD	8 %
Elongation CD	9 %
Alcohol repellent	Yes

Compliance

Clinipak RELIANCE® DEXTEX 200 BLUE/GREEN conforms to ISO 11607-1 & BS EN 868-2.

Storage/Shelf Life

It is recommended that RELIANCE® DEXTEX 200 BLUE/GREEN is stored in a cool, dry place and out of direct sunlight.

The wrap should be used within 5 years of manufacture according to the manufacture date on the product label.