



CLEAR VISON COMFORT FACE SHIELD

The Clear Vision Comfort Face Shield is designed and developed to protect the user from fluid splashes whilst providing clear vision. The longer length lens allows for enhanced user protection of the eyes, nose and mouth.

With its conformable, ergonomically designed foam headband, the Clear Vision Comfort Face Shield provides comfort and safety to users.



CLEAR VISON COMFORT FACE SHIELD – 05TD3120ST

ABOUT:

- Lens for 05TD3120ST is made of (PET) Polyethylene Terephthalate film
- Size of the film: 330mm width x 250mm length
- Forehead foam: (PU) polyurethane
- Elastic band: latex free
- Secure rivet: brass



a STERIS company

FOR MORE INFORMATION CALL 01628 810626 OR VISIT [KEYSURGICAL.CO.UK](https://www.keysurgical.co.uk)



MA-4168 Rev A, 2021