



RELIANCE HAEMOSTASIS CLIPS

RELIANCE Haemostasis Clips are used to prevent bleeding in the GI tract, which most often occurs after a polyp has been removed during a GI procedure. The clip is a small metal device that can be rotated and/or repositioned for accurate placement to join the surrounding tissue together.


ABOUT:

- Designed to clip tissue together during a GI procedure
- Repositionable and rotatable for improved accuracy, workflow, and performance
- Available in 11mm and 16mm opening widths
- Radiopaque
- MR Conditional
- Sterile
- Single-use



PRODUCT CODE	DESCRIPTION	QUANTITY
EN10400.16	RELIANCE Haemostasis Clip, Coated, 230cm length, 16mm (open width)	5/box
EN10400	RELIANCE Haemostasis Clip, Coated, 230cm length, 11mm (open width)	5/box

YOU MIGHT ALSO LIKE:



Biopsy Forceps

Polypectomy Snares



FOR MORE INFORMATION CALL 01628 810626 OR VISIT KEYSURGICAL.CO.UK

