



# AMCOR BOP

**Amcor BOP pouches and reels are a high barrier, peelable, instrument and small tray packaging solution, suitable for sterilisation.**

The range of Amcor BOP pouches and reels are made from a mix of paper and intrafilm printed film. Featuring steam sterilisation indicators and are available in flat or gusseted formats, meeting the requirements of EN 868-5 and ISO 11607-1 standards. The 60g/m2 paper offers superior resistance, high barrier and ecological options for lightweight instruments and small trays.

The BOP pouches and reels use a double barrier concept, called nesting. This system aligns the pouch or reel sizes together to allow them to slip inside one another while using the optimal amount of packaging materials.



### Sterile Compliant

The high porosity, tough packaging material and process indicators, that comply with ISO 11140-1, ensures the sterility of the device inside. The steam indicators before sterilisation are pink then once sterilised they turn to brown.



### Simple Aseptic Access

The PET and PP resistant film, decreases the risk of tearing when opening the pouch. The clean peel system minimises the release of paper fibres increasing patient safety.



### Waste Reduction

Optimised 60g/m2 paper BOP pouches and reels, minimises packaging waste and the associated cost or handling issues.



FOR MORE INFORMATION CALL 01628 810626 OR VISIT [KEYSURGICAL.CO.UK](http://KEYSURGICAL.CO.UK)





# AMCOR BOP

## Flat Pouches

### Product Number Description

90BOPS5020	50mmx200mm	90BOP22055	220mmx550mm
90BOPS8025	80mmx250mm	90BOP22060	220mmx600mm
90BOPS8015	80mmx150mm	90BOP22070	220mmx700mm
90BOPS8030	80mmx300mm	90BOP22075	220mmx750mm
90BOPS8040	80mmx400mm	90BOPS2540	250mmx400mm
90BOPS1220	120mmx200mm	90BOPS2550	250mmx500mm
90BOPS1225	120mmx250mm	90BOPS2555	250mmx550mm
90BOPS1235	120mmx350mm	90BOPS2565	250mmx650mm
90BOPS1240	120mmx400mm	90BOPS2835	285mmx350mm
90BOPS1250	120mmx500mm	90BOPS2845	285mmx450mm
90BOPS1255	120mmx550mm	90BOPS2850	285mmx500mm
90BOPS1620	160mmx200mm	90BOPS2860	285mmx600mm
90BOPS1630	160mmx300mm	90BOPS2875	285mmx750mm
90BOPS1635	160mmx350mm	90BOPS3245	320mmx450mm
90BOPS1645	160mmx450mm	90BOPS3255	320mmx550mm
90BOPS1650	160mmx500mm	90BOPS3260	320mmx600mm
90BOPS1660	160mmx600mm	90BOPS3850	380mmx500mm
90BOPS1665	160mmx650mm	90BOPS3855	380mmx550mm
90BOP22030	220mmx300mm	90BOPS3865	380mmx650mm
90BOP22040	220mmx400mm	90BOPS4260	420mmx600mm
90BOP22045	220mmx450mm	90BOPS4265	420mmx650mm

## Flat Reels

### Product Number Description

02BOPS5020	50mmx200m
02BOPS8020	80mmx200m
02BOPS1220	120mmx200m
02BOPS1620	160mmx200m
02BOPS2220	220mmx200m
02BOPS2520	250mmx200m
02BOPS2820	285mmx200m
02BOPS3220	320mmx200m
02BOPS3820	380mmx200m
02BOPS4220	420mmx200m

## Self Seal Pouches

### Product Number Description

AAPBOP7023	70mmx230mm
AAPBOP9023	90mmx230mm
AAPBOP1325	130mmx254mm
AAPBOP1336	130mmx356mm
AAPBOP2033	200mmx330mm
AAPBOP3038	300mmx390mm



FOR MORE INFORMATION CALL 01628 810626 OR VISIT [KEYSURGICAL.CO.UK](http://KEYSURGICAL.CO.UK)



Key Surgical UK is an official distributor of Amcor packaging products. Amcor is a global leader in developing and producing responsible packaging for food, beverage, pharmaceutical, medical, home and personal-care, and other products. [www.amcor.com](http://www.amcor.com) | LinkedIn