



# AMCOR ULTRA

Amcor Ultra pouches and reels are an alternative packaging solution to the traditional wrap systems used in the sterilisation process. The clear pouch is an effective and simpler way to pack medical equipment with confidence. The sterile barrier system is made out of PET film, polypropylene and polyolefin. The high strength, high porosity web provides an effective micro barrier to the packaging. The Ultra pouches are suitable for sterilisation by steam, vaporised hydrogen peroxide and ethylene oxide.



## Waste Reduction

Compared to standard non-woven wrap, the use of Ultra pouches reduces the volume of packaging materials by up to 75%.



## Time Saving

A simpler and quicker single barrier packaging system, Amcor have estimated this to be 30 seconds per tray.



## Confidence

The tough and puncture resistant pouch retains its integrity through transport and storage. This provides confidence in patient safety.



## Sterile Compliant

The higher seal strength of Ultra pouches significantly reduces the number of burst pouches during the sterilisation process. This can eliminate reprocessing and improves process efficiency.



## Simple Aseptic Access

A peel open system that eliminates the risk of tearing and does not generate fibres that could contaminate the device set. As the top film of the pouch is peeled back the device set can be easily lifted off the sterile base web of the pouch without the risk of contamination.



## Transparency

The clear film of the Ultra pouch enables users to visually confirm the contents of the device set without opening it and compromising sterility. This minimises the risk of delayed operations.



FOR MORE INFORMATION CALL 01628 810626 OR VISIT [KEYSURGICAL.CO.UK](http://KEYSURGICAL.CO.UK)






# AMCOR ULTRA

## Reel

Product Number	Description
99ULT05070	50mm x 70m
99ULT08070	80mm x 70m
99ULT12070	120mm x 70m
99ULT16070	160mm x 70m
99ULT22070	220mm x 70m
99ULT25070	250mm x 70m
99ULT28570	285mm x 70m
99ULT32070	320mm x 70m
99ULT38070	380mm x 70m
99ULT42070	420mm x 70m
99ULT47070	470mm x 70m
99ULT52070	520mm x 70m
99ULT58070	580mm x 70m


## Reel examples



## Pouches

Product Number	Description
98ULT00830	80 x 300mm
98ULT01220	120 x 200mm
98ULT01240	120 x 400mm
98ULT01630	160 x 300mm
98ULT01650	160 x 500mm
98ULT02240	220 x 400mm
98ULT02260	220 x 600mm
98ULT02550	250 x 500mm
98ULT02845	285 x 450mm
98ULT02860	285 x 600mm
98ULT03255	320 x 550mm
98ULT03845	380 x 450mm
98ULT03865	380 x 650mm
98ULT03869	380 x 690mm

## Pouches example



FOR MORE INFORMATION CALL 01628 810626 OR VISIT [KEYSURGICAL.CO.UK](http://KEYSURGICAL.CO.UK)



Key Surgical UK is an official distributor of Amcor packaging products. Amcor is a global leader in developing and producing responsible packaging for food, beverage, pharmaceutical, medical, home and personal-care, and other products. [www.amcor.com](http://www.amcor.com) | LinkedIn