




# ENDOSCOPE BRUSHES

Endoscopes are complex and often difficult to reprocess, which makes choosing the right brush critical for effective manual cleaning. Whether you need channel, valve, combination brushes, our updated line of brushes has what you need. Our single-use brushes prevent cross-contamination. Let's get to it.

## ABOUT:

Various styles available, including kits. Channel brushes feature protective tip, single use to prevent cross-contamination from re-use.



FOR MORE INFORMATION CALL 01628 810626 OR VISIT [KEYSURGICAL.CO.UK](http://KEYSURGICAL.CO.UK)





# ENDOSCOPE BRUSHES

## CHANNEL BRUSHES

PRODUCT CODE	DESCRIPTION	CHANNEL DIAM (MM)	BRUSH LENGTH	QTY
EB10106	Cleaning brush	1.2	1200mm	50/pkg
EB10136	Double-headed cleaning brush (pull through)	1.4–2.0	1200mm	100/pkg
EB10102	Double-headed cleaning brush (pull through)	2.0–4.0	2500mm	50/pkg
EB10103	Double-headed cleaning brush	2.8–4.5	2500mm	50/pkg
EN10103.3	Triple-headed cleaning brush	2.8–4.5	2500mm	50/pkg
EB10100	Double-headed cleaning brush	2.8–4.2	2500mm	50/pkg
EB10104	Double-ended cleaning brush with attached valve brush	2.8–4.2	2300mm	25/pkg
EB10105	Double-headed cleaning brush with 5mm/pipe cleaner	2.8–4.5	2300mm	50/pkg

## VALVE BRUSH

PRODUCT CODE	DESCRIPTION	CHANNEL DIAM (MM)	BRUSH LENGTH	QTY
EB10101	Biopsy valve cleaning brush	5.0 & 10	155mm	100/pkg

## ALBARRAN BRUSH

PRODUCT CODE	DESCRIPTION	CHANNEL DIAM (MM)	BRUSH LENGTH	QTY
EB10108	Designed to clean albarran channel on duodenoscope	2.3	155mm	100/pkg

## BRUSH KITS (CHANNEL BRUSHES WITH VALVE BRUSHES)

PRODUCT CODE	DESCRIPTION	CHANNEL DIAM (MM)	QTY
EN10103.3C	Triple-headed cleaning brush with valve brush (consists of EB10103.3 and EB10101)	2.8-4.5	50/pkg
EB10107C	Double-ended cleaning brush with valve brush (consists of EB10103 and EB10101)	5.0	50/pkg
EB10107G	Double-ended cleaning brush with valve brush (consists of EB10102 and EB10101)	3.0	50/pkg