

The perfect solutions to reduce wet packs and tearing...

Clinipak strives to provide solutions that influence better working practices in sterile service departments and add value to sterilisation processes.

Wet packs

The occurrences of condensation or wet packs have increased due to the growing use of high content polypropylene wraps which demonstrate advanced strength. The product selection within this pack has been designed to keep moisture near stainless steel — spreading the moisture over a larger area which allows for more efficient drying.



Protection

The reason for using of heavy duty polypropylene wraps, however, has been caused by an increase in the damage of wrapped sets, which usually occurs in the form of tears and punctures caused by sharp corners or tray feet. The additional costs required for heavy duty wraps are substantial. The products enclosed also resolve these problems in a more cost effective way.



This package contains our most popular protection product ranges which ensure maximum control.

The main benefits of these ranges include:

- Reduced costs
 - Saved essential time
 - Minimised errors
 - Lower damages
 - Fewer wet packs

PROTECTION PACK

CC001: White Rigid Card Corners

CC002: Crepe Paper Corners

TLC001: Absorbent Tray Liners

MILB022: Millennium Tray Liners

DRYMAT003: Purple Dry Mats

DRYMATGRN003: Green Dry Mats

SOK003SS: Super Soft Soaker Sheets

The standard size for each unit is 305x610mm.
Additional sizes of the samples are available on request.

Please contact our sales team for prices and enquiries:
+44 (0) 1628 810626 / sales@clinipak.co.uk

CLINIPAK