

## SCOPE BASKETS

### Product Code(s)

- CIL-xxxx
- CIL76-xxx
- SJ-102-xxx

Scope baskets are lockable stainless-steel wire mesh baskets for single endoscopes/arthroscopes with integral silicone instrument holders which are designed to protect expensive scope instruments from damage during washing and sterilising. They have an integral, lockable lid enabling the scope to be fully enclosed.



### Material

Scope Baskets are manufactured from electro-polished, 304 stainless steel with silicone instrument holders.

### Endoscope Product Codes and Sizes

➤ CIL2701	Scope Basket with small mesh Basket	640 x 150 x 77mm
➤ CIL2726	Single Scope Basket	670 x 80 x 52mm
➤ CIL2727	Single Endoscope Basket	460 x 80 x 52mm
➤ CIL2728	Single Arthroscope Basket	290 x 80 x 52mm
➤ CIL2729	Double Scope basket	460 x 80 x 52mm
➤ CIL2730	Double Scope Basket	290 x 80 x 52mm
➤ CIL2751	Double Endoscope basket	290 x 160 x 52mm
➤ CIL2753	Double Scope Basket	460 x 160 x 52mm
➤ CIL76-0001	Endoscope Basket	640 x 150 x 77mm
➤ CIL76-0010	Universal Scope Basket	580 x 160 x 52mm
➤ CIL76-0011	Universal Scope Basket	290 x 160 x 52mm
➤ CIL76-0012	Universal Scope Basket	460 x 160 x 52mm
➤ SJ-102-001P	Endoscope Basket with 1 pin	460 x 80 x 52 mm
➤ SJ-102-002P	Endoscope basket with 2 pins	460 x 80 x 52mm

### Sterilisation Compatibility

Stainless Steel Baskets are suitable for use in steam (up to 137°C) plasma sterilisation processes.