

### Product Code

- CILxxxx
- SJxxxx

A range of Stainless-Steel baskets for instrument washing, disinfection, sterilisation, storage and handling.



### Material

All baskets are made from electro-polished, 304 stainless steel (excluding CIL76-0455 which is aluminium).

### Sizes and design

➤ CIL0202	240 x 250 x 60mm	Wire Mesh Basket
➤ CIL0202LID	240 x 250mm	Lid for CIL0202 only
➤ CIL0210	480 x 250 x 60mm	Wire Mesh Basket
➤ CIL0210LID	480 x 250mm	Lid for CIL0210 only
➤ CIL0234	240 x 250 x 60mm	Flat Base Basket
➤ CIL0237	480 x 250 x 60mm	Flat Base Basket
➤ CIL0406	80 x 90 x 35mm	Fine Mesh Basket
➤ SJ/165-009	80 x 80 x 60mm	Trinket Fine Mesh basket with lid
➤ CIL0705	90 x 60 x 25mm	Trinket Box
➤ CIL0706	105 x 70 x 25mm	Trinket Box
➤ CIL0707	120 x 80 x 25mm	Trinket Box
➤ CIL24255	240 x 250 x 50mm	Flat Base Basket with Lid
➤ CIL25255	250 x 250 x 50mm	Flat Base Basket
➤ CIL25256	250 x 250 x 60mm	Flat Base Basket
➤ CIL281612	280 x 160 x 120mm	Basket with Lid
➤ CIL48256	480 x 254 x 60mm	Flat Base Basket
➤ CIL76-0030	240 x 250 x 50mm	Flat Base Basket with perforated sides
➤ CIL76-0039	380 x 240 x 50mm	Flat Base Basket with perforated sides
➤ CIL76-0047	480 x 250 x 50mm	Flat Base Basket with perforated sides
➤ CIL76-0048	480 x 250 x 60mm	Flat Base Basket with perforated sides
➤ CIL76-0049	480 x 250 x 70mm	Flat Base Basket with perforated sides
➤ CIL76-0057	480 x 360 x 60mm	Flat Base Basket
➤ CIL76-0058	480 x 480 x 50mm	Flat Base Basket
➤ CIL76-0060	480 x 480 x 70mm	Flat Base Basket
➤ CIL76-0066	520 x 250 x 50mm	Flat Base Basket with perforated sides
➤ CIL76-0076	550 x 250 x 60mm	Flat Base Basket with perforated sides
➤ CIL76-0079	550 x 430 x 50mm	Flat Base Basket with perforated sides

➤ CIL76-0086	600 x 300 x 60mm	Flat Base Basket with perforated sides
➤ SJ-131-020	650 x 250 x 70mm	Flat Base DIN basket with perforated sides
➤ CIL76-0088	670 x 300 x 70mm	Flat Base Basket
➤ CIL76-0089	700 x 250 x 60mm	Flat Base Basket
➤ CIL76-0090	700 x 300 x 90mm	Flat Base Basket with perforated sides
➤ CIL76-0095	240 x 250mm	Lid for Flat Base Basket
➤ CIL76-0096	480 x 250mm	Lid for Flat Base Basket
➤ CIL76-0143	480 x 250 x 50mm	All Wire Mesh Basket
➤ CIL76-0175	70 x 30 x 15mm	Trinket Box with perforated sides
➤ CIL76-0185	90 x 60 x 25mm	Trinket Box with wire mesh lid
➤ CIL76-0186	105 x 70 x 25mm	Trinket Box with wire mesh lid
➤ CIL76-0187	120 x 80 x 25mm	Trinket Box with wire mesh lid
➤ CIL76-0189	105 x 100 x 25mm	Trinket Box with perforated sides and base
➤ CIL76-0200	80 x 40 x 30mm	All Wire Mesh Basket
➤ CIL76-0201	80 x 80 x 35mm	All Wire Mesh Basket
➤ CIL76-0202	80 x 80 x 55mm	Fine Wire Mesh Basket
➤ CIL76-0204	105 x 70 x 25mm	Fine Wire Mesh Basket
➤ CIL76-0206	150 x 50 x 35mm	Flat Base Basket with perforated sides
➤ CIL76-0210	190 x 85 x 100mm	All Wire Mesh Basket
➤ CIL76-0215	200 x 100 x 35mm	Fine Mesh Basket
➤ CIL76-0229	220 x 140 x 50mm	Fine Mesh Basket
➤ CIL76-0256	120 x 155 x 50mm	Fine Mesh Basket, removable lid with locking clip & hinges
➤ CIL76-0270	440 x 250 x 40mm	Ophthalmic instrument basket
➤ CIL76-0475	250 x 170 x 70mm	All Wire Mesh Basket

### Sterilisation Compatibility

Clinipak economy baskets are compatible with steam sterilisation processes up to 137°C, and plasma sterilisation processes.