



ENDO-LINERS

Endo-Liners are colour-coded protective films designed to clearly identify endoscopes that are either clean or contaminated. During temporary storage or transport to and from the procedure room, we've got your scope identification needs covered - literally.

ABOUT:

- Endo-Liner (103580L) consists of two protective films that are colour-coded (red and green), one thin, opaque plastic tray liner, and one clear plastic self-adhesive tray bag.
- Green side printed with text clearly stating "CLEAN" and red side printed with text clearly stating "BIOHAZARD/CONTAMINATED" to comply with UK local requirements.
- Sterile.
- Single use.



103580L



PRODUCT CODE	DESCRIPTION	QUANTITY
103580L	Endo-Liner, Large	200/box

YOU MIGHT ALSO LIKE:

Endoscope Cleaning Brushes



FOR MORE INFORMATION CALL 01628 810626 OR VISIT KEYSURGICAL.CO.UK



MA-4155 V01