

**STAINLESS STEEL PERFORATED BASKETS**

Stainless Steel Baskets for holding instruments for sterilisation.

**Product Code(s)**

➤ BSKxxxxxxx

See attached list for available baskets sizes.

**Material**

All Baskets are manufactured using electro-polished, 304 stainless steel.

**Design**

|       |  |
|-------|--|
| Base  | Perforated Stainless-Steel Base<br>"H" Baskets have handles<br>Rounded Corners |
| Sides | Does not come with lid<br>50mm perforated Stainless-Steel sides                |

**Sterilisation Compatibility**

Stainless Steel Baskets are compatible with steam sterilisation processes up to 137°C and plasma sterilisation processes.

**SUMMARY OF STAINLESS-STEEL PERFORATED BASKET DESIGNS**

| CODE           | SIZE (mm)      | BASE TYPE  | Handles?     |
|----------------|----------------|------------|--------------|
| BSK316         | 480 X 250 X 50 | Perforated | No Handles   |
| BSK316H        | 480 X 250 X 50 | Perforated | With Handles |
| BSK317         | 250 X 250 X 50 | Perforated | No Handles   |
| BSK317H        | 250 X 250 X 50 | Perforated | With Handles |
| BSK318         | 480 X 480 X 50 | Perforated | No Handles   |
| BSK318H        | 480 X 480 X 50 | Perforated | With Handles |
| BSK580X330X75  | 580 X 330 X 70 | Perforated | No Handles   |
| BSK580X330X75M | 580 X 330 X 75 | Mesh       | No Handles   |
| BSK600X140X50  | 600 X 140 X 50 | Perforated | No Handles   |