

## CONTAINER ACCESSORIES

Containers are used worldwide to provide a protective barrier during the sterilization, handling and storage of surgical instruments. Regardless of which container model you are using, Key Surgical has just the right accessories.

Container filters are available in various shapes and sizes for a variety of container systems in accordance with the applicable standards.

Colored plastic seals in various shapes can be used for color coding and as a sustainable container closure option, paper seals are perfect for easy handling.

Container Filters


Plastic Locking Seals

Paper Locking Seals

Insert Labels

Attachment Labels

Unless otherwise stated, our products are not made with natural rubber latex.



# Container Accessories – Container Filters

## CONTAINER FILTERS

We offer a variety of filters for containers compliant with EN ISO 11607 standard. These paper filters are available with or without a process indicator and come in different sizes and shapes. Single use.








## CONTAINER FILTERS FOR AESCULAP CONTAINERS

REF	Description	Format (mm)	Qty
A	01152 round filter with steam indicator, central hole punching Ø 6mm	Ø 190	500 pcs
B	01141 round filter without indicator, central hole punching Ø 6mm, white	Ø 190	500 pcs
C	01148 small set filter with steam indicator, rectangular, with pointed corners, white	118 x 235	200 pcs
D	01186 small set filter without indicator, rectangular, with pointed corners, white	118 x 235	200 pcs
E	01317 rectangular filter with steam indicator, with pointed corners	170 x 230	500 pcs
F	01153 rectangular filter with steam indicator, with bevelled corners, white	230 x 230	500 pcs
G	01142 rectangular filter without indicator, with bevelled corners, white	230 x 230	500 pcs
H	01149 rectangular filter without indicator, with pointed corners, white	230 x 230	500 pcs

## CONTAINER FILTERS

Rectangular, white, with steam indicator. Single use.

REF	Description	Format (mm)	Qty
I	01330 rectangular filter with steam indicator, with rounded corners and a punched hole (long)	104 x 104	250 pcs
J	01158 rectangular filter with steam indicator, with sharp corners and two punched holes (round)	110 x 240	200 pcs
K	01143 rectangular filter with steam indicator, with pointed corners	180 x 270	500 pcs
L	01321 rectangular filter with steam indicator, with pointed corners	180 x 255	500 pcs






Container Accessories

## Container Accessories – Container Filters


### CONTAINER FILTERS

Rectangular, white. Compliant with EN ISO 11607 standard. Single use.

REF	Description	Format (mm)	Qty
A 01144	without indicator, rectangular, pointed corners	240 x 240	250 pcs
A 01198	without indicator, rectangular, pointed corners	230 x 540	250 pcs
B 01304	without indicator, rectangular, 4 punched holes, with rounded corners	210 x 490	100 pcs
B 01318	without indicator, rectangular, 4 punched holes, with rounded corners	132 x 358	50 pcs
C 01346	with steam indicator	300 x 300	500 pcs
C 01347	with steam indicator	300 x 600	500 pcs




## Container Accessories – Container Filters



A



B



C

### CONTAINER FILTERS

Non-woven container filters, 67gsm according to EN ISO 11607. Round or square format. Single use.

REF	Description	Format (mm)	Qty
A 01145	square, green, with steam indicator, pointed corners, 2 punched holes	225 x 225	500 pcs
B 01378	round, blue, for H <sub>2</sub> O <sub>2</sub> sterilization, with steam indicator	190 x 6.0	500 pcs
01365	round, blue, with steam indicator	190 x 6.0	500 pcs
C 01368	round, green, without indicator, filter imprint: blank	Ø 92	1000 pcs
01188	pointed corners, blue, for H <sub>2</sub> O <sub>2</sub> sterilization, without indicator	120 x 235	500 pcs
01189	pointed corners, blue, for H <sub>2</sub> O <sub>2</sub> sterilization, without indicator	130 x 110	500 pcs



## Container Accessories – Container Filters

### CONTAINER FILTERS FOR OTHER CONTAINERS

Square and rectangular, white, with indicator. Single use.

REF	Description	Format (mm)	Qty
A	01314 rectangular filter, with indicator, with pointed corners	145 x 205	100 pcs
B	01315 rectangular filter, with indicator, with bevelled corners	240 x 240	250 pcs
C	01327 rectangular filter, with indicator, with pointed corners	240 x 540	100 pcs

### CONTAINER FILTERS FOR ALL WAGNER CONTAINERS


Round and rectangular, white, with indicator. Single use.

REF	Description	Format (mm)	Qty
F	01355 round, central hole punching Ø 10mm	Ø 190	500 pcs
G	01369 round, central hole punching Ø 10mm	Ø 75	1000 pcs
H	01363 rectangular, with pointed corners	245 x 400	100 pcs
I	01323 rectangular, with pointed corners	96 x 212	250 pcs

### STAINLESS STEEL DISPENSER FOR FILTERS

Including adhesive material for wall mounting.

REF	Format (mm)	Qty
J	20437 195 x 195 x 100	1 pcs




A



B




F



C



G



H







I



J

## Container Accessories – Plastic Locking Seals



### PLASTIC LOCKING SEALS

We offer a wide range of different designs to suit all type of containers. Our seals are available both with and without process indicators. They are used as a tamper free lock and have a predetermined breaking point for easy opening. We have options for all common container makes and a wide range of colors. We also offer numbered security seals (REF 02639) for the sealed return of loan instruments, implants, and containers. Single use.



### LOCKING SEAL - ANGULAR BRIDGE

Bridge length 19mm.

REF	Description	Color	Qty	
A	01020	without indicator	dark blue	1000 pcs
B	01021	without indicator	green	1000 pcs
C	01022	without indicator	white	1000 pcs
D	01023	without indicator	red	1000 pcs
E	01024	without indicator	yellow	1000 pcs
F	01234	without indicator	orange	1000 pcs
H	01286	without indicator	black	1000 pcs
I	02623	without indicator	light blue	1000 pcs
J	01228	with steam indicator	dark blue	1000 pcs
K	01229	with steam indicator	green	1000 pcs
L	01230	with steam indicator	white	1000 pcs
M	01231	with steam indicator	red	1000 pcs
N	01232	with steam indicator	yellow	1000 pcs
O	01235	with steam indicator	orange	1000 pcs
P	01294	with steam indicator	black	1000 pcs
Q	02670	with steam indicator	light blue	1000 pcs



## Container Accessories – Plastic Locking Seals



### UNIVERSAL LOCKING SEALS UNILOCK - ANGULAR BRIDGE

Bridge length 19mm. Single use.

REF	Description	Color	Qty
A 01240	without indicator	blue	1000 pcs
B 01241	without indicator	green	1000 pcs
C 01242	without indicator	white	1000 pcs
D 01244	without indicator	yellow	1000 pcs
E 01245	with steam indicator	blue	1000 pcs
F 01246	with steam indicator	green	1000 pcs
G 01247	with steam indicator	white	1000 pcs
H 01248	with steam indicator	red	1000 pcs
I 01249	with steam indicator	yellow	1000 pcs



### UNIVERSAL LOCKING SEALS UNILOCK XL - SQUARE BRIDGE

Bridge length 27mm. Single use.

REF	Description	Color	Qty
J 01250	without indicator	blue	1000 pcs
K 01251	without indicator	green	1000 pcs
L 01252	without indicator	white	1000 pcs
M 01254	without indicator	yellow	1000 pcs
N 02514	without indicator	violet	1000 pcs
O 02515	without indicator	orange	1000 pcs
P 01255	with steam indicator	blue	1000 pcs
Q 01256	with steam indicator	green	1000 pcs
R 01257	with steam indicator	white	1000 pcs
S 01258	with steam indicator	red	1000 pcs
T 01259	with steam indicator	yellow	1000 pcs



## Container Accessories – Plastic Locking Seals



### UNIVERSAL LOCKING SEALS - ANGULAR BRIDGE

Suitable for Aesculap and Wagner containers. Bridge length 34mm. Single use.

REF	Description	Color	Qty
A 01025	without indicator	dark blue	1000 pcs
B 01026	without indicator	red	1000 pcs
C 01027	without indicator	green	1000 pcs
D 01028	without indicator	yellow	1000 pcs
E 01029	without indicator	white	1000 pcs
F 02624	without indicator	orange	1000 pcs
G 02625	without indicator	violet	1000 pcs
H 01220	with steam indicator	dark blue	1000 pcs
I 01221	with steam indicator	red	1000 pcs
J 01233	with steam indicator	green	1000 pcs
K 01238	with steam indicator	yellow	1000 pcs
L 01239	with steam indicator	white	1000 pcs

### 'AW' SEALS - ROUND BAR


Suitable for Aesculap and Wagner containers. Bridge length 19mm. Single use.



REF	Description	Color	Qty
M 02612	without indicator	red	1000 pcs
N 02613	without indicator	blue	1000 pcs
O 02614	without indicator	green	1000 pcs
P 02615	without indicator	yellow	1000 pcs
Q 02661	without indicator	white	1000 pcs
R 02662	without indicator	black	1000 pcs
S 02632	with steam indicator	red	1000 pcs
T 02633	with steam indicator	blue	1000 pcs
U 02634	with steam indicator	green	1000 pcs
V 02635	with steam indicator	yellow	1000 pcs



## Container Accessories – Plastic Locking Seals



### UNIVERSAL SEALS - ROUND BRIDGE

Length 195mm. Single Use.


REF	Description	Color	Qty
<b>A</b> 01204	without indicator	dark blue	1000 pcs
<b>B</b> 02643	without indicator	yellow	1000 pcs
<b>C</b> 01205	with steam indicator	dark blue	1000 pcs

### UNIVERSAL SEALS - ANGULAR BRIDGE

Bridge length 80mm. Single use.

REF	Description	Color	Qty
<b>D</b> 01014	without indicator	blue	1000 pcs
<b>E</b> 01216	without indicator	red	1000 pcs
<b>F</b> 01217	without indicator	green	1000 pcs
<b>G</b> 01218	without indicator	yellow	1000 pcs
<b>H</b> 01219	without indicator	black	1000 pcs
<b>I</b> 01013	with steam indicator	blue	1000 pcs
<b>J</b> 01212	with steam indicator	red	1000 pcs
<b>K</b> 01213	with steam indicator	green	1000 pcs
<b>L</b> 01214	with steam indicator	yellow	1000 pcs
<b>M</b> 01215	with steam indicator	black	1000 pcs

## Container Accessories – Plastic Locking Seals



### PLASTIC SEALS FOR MARTIN CONTAINERS

REF	Description	Color	Qty
A 01011	without indicator	white	1000 pcs
B 01012	with steam indicator	white	1000 pcs

### SECURITY SEALS

REF	Description	Color	Qty
C 02639	290mm long, usable length 220mm	yellow	400 pcs



## Container Accessories – Printed Locking Tags



### LOCKING TAGS

Versatile, durable, and economical. Tamper-resistant. Compatible with steam and H<sub>2</sub>O<sub>2</sub> sterilization. Single use. Customizable. 25.4mm length.



REF	Details	Color	Qty
A LT-100N*	uniquely numbered	white	100 pcs
LT-100N-1000*			1000 pcs
B LT-300N*	uniquely numbered	green	100 pcs
LT-300N-1000*			1000 pcs
C LT-400N*	uniquely numbered	red	100 pcs
LT-400N-1000*			1000 pcs
D LT-500N*	uniquely numbered	yellow	100 pcs
LT-500N-1000*			1000 pcs
E LT-900N*	uniquely numbered	black	100 pcs
LT-900N-1000*			1000 pcs






## PAPER SEALS

We offer you an environmentally friendly way of closing containers. The seals can be attached to almost all common containers and follow a simple principle. With our paper seals you ensure safe closure for transport and storage due to thermal closure of the components. Available in different versions with and without documentation fields. In some cases, it is also possible to stick on labels. Single use.

REF	Description	Color	Format (mm)	Qty
A	01010	white	100 x 18	1000 pcs
B	01003	white	105 x 30	1000 pcs
C	01004	white	117 x 38	500 pcs
D	01005	white	117 x 38	500 pcs
E	02638	white	105 x 30	1000 pcs
F	04120	white	132 x 38	500 pcs
G	04179	white	145 x 35	750 pcs
H	01006	white	35 x 74	1000 pcs
I	04114	white	35 x 74	1000 pcs
J	01008	white	35 x 74	1000 pcs
K	01009	white	35 x 74	1000 pcs
L	04108	white	35 x 74	1000 pcs
M	04113	white	35 x 80	1000 pcs
61006	closure label, with indicator, self-adhesive, insert from the right, with documentation field	white	74 x 35	1000 pcs






A

B

C



D

E

F



G

H

**INSERT LABELS**

Insert labels make documentation easier for you, they can simply be pushed into the opening provided on the container. Available in different designs suitable for Aesculap, Martin, or other container systems. We also offer blank versions and with documentation fields in which sterilization data, batch numbers, and sieve names can be noted. Our insert labels are made of special cardboard compatible with sterilization. Single use.

**DOCUMENTATION LABELS FOR WAGNER CONTAINERS**


With indicator and documentation field.

REF	Self-adhesive	Format (mm)	Qty
A 01133	no	40 x 59	1000 pcs
B 01134	yes	40/10 x 59	1000 pcs
C 01169	yes	40 x 59/74	1000 pcs

**DOCUMENTATION LABELS FOR MISCELLANEOUS CONTAINERS**

REF	Description	Format (mm)	Qty
D 04067	with indicator, with documentation field, self-adhesive; for Geuder containers	23 x 60	1000 pcs
E 04116	with indicator, with documentation field, not self-adhesive; for CBM containers	39 x 55	1000 pcs
F 04112	with indicator, with documentation field, self-adhesive; for Microstop containers	39 x 63	1000 pcs
G 01116	with indicator, with documentation field, not self-adhesive; for sterilization containers	38 x 67	1000 pcs
H 04003	set label with indicator, self-adhesive, carrier material slotted from special cardboard 280gsm	40 x 80	1000 pcs
I 04105	without indicator, blank, self-adhesive; for sterilization containers	46 x 60	1000 pcs





**DOCUMENTATION LABELS FOR AESCULAP CONTAINERS**

For single use.

REF	Description	Format (mm)	Qty
A	04102 blank, not self-adhesive	39 x 68	1000 pcs
B	01400 with indicator, with documentation field, not self-adhesive	35 x 80	1000 pcs
C	04184 without indicator, blank, not self-adhesive	35 x 80	1000 pcs
D	01135 with indicator, with documentation field, not self-adhesive	39 x 68	1000 pcs
E	01136 with indicator, with documentation field, self-adhesive	39 x 67	1000 pcs
F	04000 with indicator, with documentation field, self-adhesive	18 x 60	1000 pcs
G	04005 with indicator, with documentation field, not self-adhesive, 100 pieces each packed	18 x 50	1000 pcs
	04107 with indicator 'STERILE', without frame, self-adhesive	39 x 67	1000 pcs

**DOCUMENTATION LABELS FOR MARTIN CONTAINERS**

With steam indicator and documentation field. Single use.

REF	Self-adhesive	Format (mm)	Qty
H	yes	46 x 50	1000 pcs
I	yes	46 x 60	1000 pcs

**DOCUMENTATION LABELS FOR MARTIN CONTAINERS**

For single use.

REF	Description	Format (mm)	Qty
J	01120 with indicator, with documentation field, self-adhesive	36 x 67	1000 pcs
K	04123 blank, not self-adhesive	36 x 67	1000 pcs
L	01121 with indicator, with documentation field, not self-adhesive	36 x 64	1000 pcs



## Container Accessories – Tag Labels


### TAG LABELS

Key Surgical tags help you to quickly label your containers and trolleys. They are heat-resistant up to 134 °C, not self-adhesive. Single use.


REF	Description	Format (mm)	Qty
A	04162 with hole, white, made of polyester	90 x 50 x 20	1000 pcs
B	04192 with hole, red on one side, made of polyester	90 x 50 x 20	1000 pcs
C	04193 with hole, blue on one side, made of polyester	90 x 50 x 20	1000 pcs
D	04194 with hole, yellow on one side, made of polyester	90 x 50 x 20	1000 pcs
E	04195 with hole, red on both sides, made of polyester	90 x 50 x 20	1000 pcs
	04216 with hole, white, made of polyester	110 x 60 x 20	1000 pcs
	04289 with hole, red on both sides, made of polyester	110 x 60 x 20	1000 pcs
G	04172 with cross-cut, white, made of polyester	75 x 40 x 20	1000 pcs



A




B



C



D



E



G