

CSSD EQUIPMENT

On the following pages you will find various products for equipping your CSSD.

The importance of hospital logistics should not be underestimated when it comes to the reprocessing of medical devices and instruments. Short and long distances, safe transport and the avoidance of nosocomial (hospital-acquired) infections are all issues which need to be taken into account.

Dustproof bags protect the sterile barrier, making them perfectly suited for transporting and storing sterile supplies. Key Surgical provides dustproof bags in various sizes and material strengths and with various sealing systems.

Our packing table sheets made of silicone serve, among other things, as impact protection for instruments, and they also make handling easier thanks to the non-slip properties and they also lessen the unpleasant cold, glare, and noise effects of metal surfaces by staff of your CSSD.

[Transport and Disposal](#)

[Dustproof Bags](#)

[Rinsing Guns](#)

[Packing Table Sheets](#)

Unless otherwise stated, our products are not made with natural rubber latex.



A DISPOSAL BOXES MADE OF POLYPROPYLENE

Uncolored with clips for color coding, suitable for washer-disinfectors and leak-proof.

REF	Description	Lid	Drain	Format (mm)	Capacity (litre)	Qty
09811	for washer-disinfectors	transparent	no	325 x 256 x 100	4	1 pc
09813	for washer-disinfectors	transparent	no	530 x 325 x 100	8	1 pc
09824	for washer-disinfectors	without	no	530 x 325 x 150	15	1 pc
09844	for washer-disinfectors	transparent	no	530 x 325 x 200	20	1 pc

A



DISPOSAL BOXES MADE OF POLYPROPYLENE

Liquid-tight, with drip tray insert, lid can be sealed on all sides, suitable for washer-disinfectors.

REF	Lid	Drop-sieve	Format (mm)	Capacity (litre)	Qty
A 09816	transparent	yes	278 x 115 x 103	2	1 pc
A 09817	transparent	yes	297 x 152 x 118	3.4	1 pc
A 09806	transparent	yes	295 x 230 x 106	4.8	1 pc
A 09807	transparent	yes	295 x 230 x 138	6.5	1 pc
A 09818	transparent	yes	295 x 230 x 203	10	1 pc
A 09890	transparent	no	361 x 274 x 212	16	1 pc
A 09891	transparent	no	456 x 300 x 225	21	1 pc



A



DISPOSAL AND DISINFECTION TRAYS AND ACCESSORIES MADE OF POLYPROPYLENE

These containers are used for the safe transport of surgical instruments and also for longer transport routes, for example when reprocessing outside the facility. With reusable seals, the tub can be securely closed with the lid. Various cable ties and security seals can be attached for color coding. Compatible for washer-disinfectors, not autoclavable, heat-resistant up to 110 °C.

REF	Description	Color	Suitable for REF	Format (mm)	Qty
A 17209	disposal tray	yellow	17208 17210	590 x 295 x 150	1 pc
A 17211	disposal tray	white	17208 17210	590 x 295 x 150	1 pc
A 17208	lid for disposal and disinfection tub	yellow	17209 17211	590 x 295 x 30	1 pc
A 17210	lid for disposal and disinfection tub	white	17209 17211	590 x 295 x 30	1 pc
A 17212	security seal for disposal tray with lid	white	17208 17209 17210 17211	-	10 pcs
A 17213	security seal for disposal tray with lid	white	17208 17209 17210 17211	-	100 pcs

B



SOIL SIEVES MADE OF STAINLESS STEEL

Ideal as drip inserts for trays.

REF	Description	Format (mm)	Qty
B 08956	with handles, mesh size 4.2 x 4.2	245 x 175 x 90	1 pc

POLYPROPYLENE CONTAINERS

Transparent, made of polypropylene. Security seals available. Heat-resistant up to 75 °C, not autoclavable.

REF	Lid	Outer dimensions (mm)	Inner dimensions (mm)	Qty
A BGDNN71100B	blue	320 x 260 x 100	285 x 252 x 100	1 pc
BGDNN71100G	yellow	320 x 260 x 100	285 x 252 x 100	1 pc
BGDNN71100R	red	320 x 260 x 100	285 x 252 x 100	1 pc
BGDNN70100B	blue	525 x 325 x 100	483 x 285 x 100	1 pc
BGDNN70100G	yellow	525 x 325 x 100	483 x 285 x 100	1 pc
BGDNN70100R	red	525 x 325 x 100	483 x 285 x 100	1 pc
BGDNN70100	transparent	525 x 325 x 100	483 x 285 x 100	1 pc
BGDNN70200G	yellow	525 x 325 x 200	483 x 285 x 200	1 pc
BGDNN70200R	red	525 x 325 x 200	483 x 285 x 200	1 pc
BGDNN70200	transparent	525 x 325 x 200	483 x 285 x 200	1 pc



A



B



CART IDENTIFICATION BOARD

Suitable for washer-disinfectors and with rounded corners.

REF	Description	Front	Back	Format (mm)	Qty
17996	-	green with imprint: 'aufbereitete Medizinprodukte'	red with imprint: 'GEBRAUCHTES MEDIZINISCHES INSTRUMENT'	210 x 297	10 pcs
17997	with information on UN regulation	green with imprint: 'aufbereitete Medizinprodukte'	red with imprint: 'GEBRAUCHTES MEDIZINISCHES INSTRUMENT'	210 x 297	10 pcs
A 17998	-	green with imprint: 'aufbereitete Medizinprodukte'	red with imprint: 'GEBRAUCHTES MEDIZINISCHES INSTRUMENT'	210 x 297	10 pcs
B 17999	with information on UN regulation	green with imprint: 'aufbereitete Medizinprodukte'	red with imprint: 'USED MEDICAL INSTRUMENT'	210 x 297	10 pcs



DUSTPROOF BAGS WITH FLAP AND ADHESIVE SEAL

Solid design LDPE flat pouches. Transparent. Single use, not suitable for sterilization.

REF	Description	Format (mm)	Material thickness μ	Qty
A 03004	unprinted	100 x 150	50	500 pcs
A 03023	unprinted	150 x 270	50	100 pcs
A 03025	unprinted	150 x 380	50	100 pcs
A 03029	unprinted	200 x 300	50	100 pcs
A 03028	unprinted	200 x 650	80	100 pcs
A 03007	unprinted	300 x 400	60	100 pcs
A 03030	unprinted	300 x 500	50	100 pcs
A 03032	unprinted	600 x 750	50	100 pcs
A 03037	unprinted	630 x 820	60	100 pcs
A 03049	unprinted	700 x 1050	50	50 pcs
B 03033	unprinted	245 x 345	60	100 pcs
B 03036	unprinted	300 x 600	60	100 pcs
B 03017	unprinted	400 x 600	50	100 pcs
B 03031	unprinted	420 x 700	50	100 pcs
B 03022	unprinted	450 x 480	30	100 pcs
B 03035	unprinted	460 x 750	60	100 pcs
B 03034	unprinted	640 x 770	60	100 pcs
B 03040	with printing in German, English, Italian and French	500 x 500	60	100 pcs

DUSTPROOF BAGS WITH SLIDER, RESEALABLE - UNPRINTED

Solid design made of LDPE, flat pouch transparent. Single use. Not suitable for sterilization.

REF	Format (mm)	Material thickness μ	Qty
C 03041	320 x 230	75	100 pcs
C 03046	480 x 640	70	50 pcs
C 03045	320 x 440	70	100 pcs
C 03044	170 x 120	60	100 pcs
C 03043	200 x 150	60	100 pcs

MICROSAP BAGS - RESEALABLE - UNPRINTED

Solid design made of LDPE, flat pouch. Material thickness 50 µ. Transparent. Single use, not suitable for sterilization.



REF	Format (mm)	Qty
A 09919	90 x 670	1000 pcs
A 09906	100 x 150	1000 pcs
A 09913	120 x 170	1000 pcs
A 09903	150 x 220	1000 pcs
A 09914	180 x 250	1000 pcs
A 09900	220 x 310	1000 pcs
A 09905	300 x 400	1000 pcs
A 09911	400 x 500	500 pcs



A

CSSD Equipment – Dustproof Bags

LDPE SECURITY BAG

For the hygienic storage after disinfection. The security closure is 40mm high and the tear-off label 50mm, with serial number with anti-tamper closure. For single use.

REF	Format (mm)	Internal dimensions (mm)	Qty
A 03077	520 x 500	500 x 500	100 pcs
A 03078	770 x 500	750 x 500	100 pcs

FLAT POUCHES / DUSTPROOF BAGS - UNPRINTED

For single use.

REF	Description	Format (mm)	Material thickness μ	Qty
B 09916	cleanroom quality, in secondary packaging	70 x 700	50	250 pcs
C 09924	brown-opaque color	200 x 400	80	1000 pcs



B



C



RINSING GUN

The stainless steel cleaning gun is sterilizable and easy to disassemble. The cleaning gun was developed to clean, rinse and dry off instruments. The Hela Cleaning Kit is ideal for cleaning and drying flexible endoscopes. The cleaning gun may be used with demineralized water (deionized water), tap water or compressed air.

CLEANING GUNS FOR HELA

	REF	Model	Hose	Connection	Qty
A	49994	Hela	1.5 m stainless steel hose	1/2"	1 pc
A	49992	Hela	3.0 m stainless steel hose	1/2"	1 pc
A	49993	Hela	1.5 m stainless steel hose	3/4"	1 pc
A	49991	Hela	3.0 m stainless steel hose	3/4"	1 pc
B	49998	Hela	1.5 m PVC fabric hose	1/2"	1 pc
B	49996	Hela	3.0 m PVC fabric hose	1/2"	1 pc
B	49997	Hela	1.5 m PVC fabric hose	3/4"	1 pc
B	49995	Hela	3.0 m PVC fabric hose	3/4"	1 pc

C HELA CLEANING-KIT WITH PREMIUM SET

Cleaning kit does not include a cleaning gun.


	REF	Model	Set consists of:	Hose	Connection	Qty
	49984	Hela	6 screwable attachments angle holder	stainless steel	1/2"	1 pc
	49986	Hela	6 screwable attachments angle holder	stainless steel	1/2"	1 pc
	49985	Hela	6 screwable attachments angle holder	stainless steel	3/4"	1 pc
	49990	Hela	6 screwable attachments angle holder	1.5 m PVC fabric	1/2"	1 pc
	49983	Hela	6 screwable attachments angle holder	stainless steel	3/4"	1 pc
	49988	Hela	6 screwable attachments angle holder	3.0 m PVC fabric	1/2"	1 pc
	49989	Hela	6 screwable attachments angle holder	1.5 m PVC fabric	3/4"	1 pc
	49987	Hela	6 screwable attachments angle holder	3.0 m PVC fabric	3/4"	1 pc
	49982	Hela	6 screwable attachments angle holder	without	without	1 pc



RINSING GUNS

REF	Model	Description	Connection	Qty
A 49999	Dina	can be disassembled, made of aluminum	without	1 pc
A 49975	Dina	1.5 m stainless steel hose	1/2"	1 pc
A 49973	Dina	3 m stainless steel hose	1/2"	1 pc
A 49974	Dina	1.5 m stainless steel hose	3/4"	1 pc
A 49972	Dina	3 m stainless steel hose	3/4"	1 pc
B 49979	Dina	1.5 m PVC fabric hose	1/2"	1 pc
B 49977	Dina	3 m PVC fabric hose	1/2"	1 pc
B 49978	Dina	1.5 m PVC fabric hose	3/4"	1 pc
B 49976	Dina	3 m PVC fabric hose	3/4"	1 pc
C 49963	Nova	rigid, made of aluminum	without	1 pc
C 49958	Nova	1.5 m stainless steel hose	1/2"	1 pc
C 49956	Nova	3 m stainless steel hose	1/2"	1 pc
C 49957	Nova	1.5 m stainless steel hose	3/4"	1 pc
C 49955	Nova	3 m stainless steel hose	3/4"	1 pc
D 49962	Nova	1.5 m PVC fabric hose	1/2"	1 pc
D 49960	Nova	3 m PVC fabric hose	1/2"	1 pc
D 49961	Nova	1.5 m PVC fabric hose	3/4"	1 pc
D 49959	Nova	3 m PVC fabric hose	3/4"	1 pc

A



A



A CLEANING SET WITH HOSES

REF	Model	Set consists of:	Hose	Connection	Qty
49967	Dina	8 attachments and holders made of stainless steel	stainless steel	1/2"	1 pc
49965	Dina	8 attachments and holders made of stainless steel	stainless steel	1/2"	1 pc
49966	Dina	8 attachments and holders made of stainless steel	1.5 m stainless steel hose	3/4"	1 pc
49964	Dina	8 attachments and holders made of stainless steel	3.0 m stainless steel hose	3/4"	1 pc
49971	Dina	8 attachments and holders made of stainless steel	1.5 m PVC fabric	1/2"	1 pc
49969	Dina	8 attachments and holders made of stainless steel	3 m PVC fabric	1/2"	1 pc
49970	Dina	8 attachments and holders made of stainless steel	1.5 m PVC fabric	3/4"	1 pc
49968	Dina	8 attachments and holders made of stainless steel	3 m PVC fabric	3/4"	1 pc
49950	Nova	8 attachments and holders made of stainless steel	1.5 m stainless steel hose	1/2"	1 pc
49948	Nova	8 attachments and holders made of stainless steel	3 m stainless steel hose	1/2"	1 pc
49949	Nova	8 attachments and holders made of stainless steel	1.5 m stainless steel hose	3/4"	1 pc
49947	Nova	8 attachments and holders made of stainless steel	3 m stainless steel hose	3/4"	1 pc
49954	Nova	8 attachments and holders made of stainless steel	1.5 m PVC fabric	1/2"	1 pc
49952	Nova	8 attachments and holders made of stainless steel	3 m PVC fabric	1/2"	1 pc
49953	Nova	8 attachments and holders made of stainless steel	1.5 m PVC fabric	3/4"	1 pc
49951	Nova	8 attachments and holders made of stainless steel	3 m PVC fabric	3/4"	1 pc




CSSD Equipment – Packing Table Sheets

PACKING TABLE SHEETS

Reprocessing cycles are getting shorter and cycle rates are often increasing due to reductions in staffing levels. With these pressures, inattentiveness when handling instruments is a real risk. Our packing table sheets protect instruments and also aid handling thanks to their non-slip properties. Moreover, they reduce the unpleasant cold and glare of metallic surfaces for CSSD staff. The sheets are available in all formats up to a maximum of 300cm (length) x 120cm (width) and come in a range of colors. The sheets are 2mm thick and are both highly robust and durable. Suitable for washer-disinfectors, heat-resistant.

REF	Color	Thickness (mm)	Qty
A 08700	transparent	2	1 pc
B 08705	red-brown	2	1 pc
C 08703	blue	2	1 pc
D 08701	leaf green	2	1 pc

Transport and Disposal – Transportation Boxes



TRANSPORTATION BOXES

Stackable boxes that can be custom printed with text on the box and lid. They can withstand washer-disinfectors for decontamination purposes. BTB005 must be used with a trolley BTB005TROL.

- BTB001 suitable for transporting packing group I dangerous substances
- BTB002 suitable for transporting packing group II dangerous substances
- BTB003 suitable for transporting packing group I dangerous substances
- BTB004A suitable for transporting packing group II dangerous substances
- BTB005 suitable for transporting packing group II dangerous substances

REF	Weight (Kg)	Capacity (liters)	Color	Size (external)	Qty
A BTB001	0.6	4	blue	300 x 200 x 120mm	1 pc
B BTB002	0.9	7	blue	300 x 200 x 170mm	1 pc
C BTB003	1.4	16	blue	400 x 300 x 220mm	1 pc
D BTB004A	3.3	49	blue	600 x 400 x 320mm	1 pc
BTB005	10.2	170	blue	800 x 600 x 600mm	1 pc
E RTB002	0.9	7	red	300 x 200 x 170mm	1 pc
F RTB003	1.4	16	red	400 x 300 x 220mm	1 pc
G RTB004A	3.3	49	red	600 x 400 x 320mm	1 pc

INNER BASKETS

Powder-coated aluminum inner baskets designed to fit transportation boxes.

REF	Description	Qty
H BTB002BASK	inner basket to fit BTB002 and RTB002	1 pc
BTB003BASK	inner basket to fit BTB003 and RTB003	1 pc
BTB004BASK	inner basket to fit BTB004A and RTB004A	1 pc

Transport and Disposal – Transportation Boxes



TRANSPORTATION BOX OPENER

Stainless steel opener suitable for all BTB and RTB boxes.

REF	Description	Qty
A	BTBOPENER box opener for BTB and RTB boxes, stainless steel	1 pc

TRANSPORTATION BOX TROLLEYS

Painted steel trolley with wheels for BTB and RTB boxes. Suitable for washer-disinfectors.

REF	Description	Qty
B	BTB004TROL trolley to fit BTB004A and RTB004A	1 pc
C	BTB005TROL trolley to fit BTB005	1 pc




TRANSPORTATION BOX LABELS

A selection of labels for 'BIOHAZARD' or 'USED MEDICAL DEVICE' to be inserted in plastic sleeves of transportation boxes. Our biohazard labels conform to the requirements for a 6.2 hazard infectious substances.

REF	Description	Size (mm)	Qty
D	BTBLABEL001 biohazard, non-laminated label	125 x 40	20 pcs
	BTBLABEL003 biohazard, laminated label	125 x 40	20 pcs
	BTBLABEL004 biohazard, laminated label	150 x 70	20 pcs
E	BTB232SML used medical device, laminated label	125 x 40	100 pcs
E	BTB232LRG used medical device, laminated label	150 x 70	100 pcs





SECURITY CABLE TIES

Plastic cable ties in a range of colors for use during transportation and identification.

REF	Color	Size (mm)	Qty
A BTBREDTIE	red	3.6 x 145	100 pcs
B BTBYELTIE	yellow	3.6 x 145	100 pcs
C TIE001	blue	3.5 x 145	100 pcs
D TIE006	clear	1.5 x 125	500 pcs
E TIE010	red	4.8 x 300	100 pcs

CABLE TIES

Plastic cable ties in a range of colors for use during transportation. Colored cable ties help identification, along with numbered cable ties.

REF	Description	Color	Size (mm)	Qty
F TIE003	individually numbered	green	2 x 200	100 pcs
TIE005	individually numbered	red	2 x 200	100 pcs