

DOCUMENTATION

Ever stricter standards mean that complete documentation is now part of everyday routine when reprocessing instruments. Key Surgical offers a vast array of labels for identification of your sterile goods. Our range includes both paper labels for single use and longlasting foil-based labels. Take advantage of our comprehensive range and label your medical devices to meet your specific requirements.

Our foil and paper labels are available in almost all formats, and for paper labels you can choose between those with and without indicators for reliable identification in ongoing production for all common sterilization processes.

In our product range you will also find labels for low-temperature processes, such as H₂O₂ and thermo transfer foil I in a selection of qualities.

Paper Labels

Foil Labels

Thermo Transfer Foil

Pre-printed Labels Labellers and Labelling Accessories

Miscellaneous Documentation Items

Unless otherwise stated, our products are not made with natural rubber latex.



Document the contents of your packaged medical devices using our high-quality paper labels. The labels are compatible with virtually all label printers, so you can print them according to your requirements and are thus always sure of finding the most important information at a glance. In the Key Surgical portfolio you will find labels in common formats. The duplex self-adhesive labels are particularly suitable for documentation in patient files. You will also find in our portfolio labels wound on the inside or outside and on various reel cores. Labels without tear-off perforations are ideal for cut-off printers.

SINGLE LABELS AND ENDLESS LABELS - WITHOUT INDICATOR

Blank without indicator, wound outside. For single use.

REF	Description	Self-adhesive	Tear-off perforation	Format (mm)	Core-Ø (mm)	Qty
A	roll of 1000 labels	duplex	yes	50 x 33	40	1 roll
B	suitable for non-woven materials, roll of 1000 labels	duplex	yes	70 x 33	40	1 roll
C	roll of 1000 labels	duplex	yes	60 x 34	40	1 roll
D	divided into 3 individual labels, roll of 2000 labels	duplex	yes	102 x 34 51 x 22 51 x 22	76	1 roll
D	divided into 3 individual labels, roll of 500 labels	duplex	yes	102 x 34 51 x 22 51 x 22	40	1 roll
E	roll of 2000 labels	duplex	yes	75 x 35	76	1 roll
E	with slots, heat-resistant up to 134 °C, roll of 2000 labels	duplex	no	75 x 35	76	1 roll
E	roll of 2000 labels	duplex	no	63 x 40	76	1 roll



A



B




C



D



E




A



B




C




D



E



F




G


SINGLE LABELS - WITH STEAM INDICATOR (PINK TO BROWN)

Wound on the outside and duplex self-adhesive. Blank with steam indicator.

REF	Description	Tear-off perforation	Format (mm)	Core-Ø (mm)	Qty
04322	upper material with horizontal slits - roll of 1000 labels	yes	60 x 33	25	1 roll
A 04010	with red colored corner - roll of 1000 labels	yes	60 x 33	40	1 roll
B 04086	roll of 1000 labels	yes	60 x 33	40	1 roll
B 04158	roll of 2000 labels	yes	60 x 33	40	1 roll
C 04214	roll of 1000 labels	yes	59 x 35	40	1 roll
D 04209	upper material with horizontal slits - roll of 1000 labels	yes	59 x 35	40	1 roll
E 04290	roll of 1000 labels	yes	58 x 38	40	1 roll
F 04141	roll of 1000 labels	yes	57 x 40	40	1 roll
G 04167	with slit in the upper material - roll of 660 labels	yes	72 x 48	40	1 roll
04287	with slit in the upper material - roll of 660 labels	yes	72 x 48	76	1 roll



A




B




C




D



E



G



H

SINGLE LABELS - WITH STEAM INDICATOR SUITABLE FOR NON-WOVEN MATERIAL

Blank with steam indicator (pink to brown). Wound outside and duplex self-adhesive.


REF	Description	Tear-off perforation	Format (mm)	Core-Ø (mm)	Qty
A	04160 roll of 1000 labels	yes	60 x 31	25	1 roll
B	04207 roll of 1000 labels	yes	50 x 33	40	1 roll
C	04089 roll of 1000 labels	no	70 x 33	40	1 roll
D	04301 with information on standards and color change roll of 1000 labels	yes	70 x 33	40	1 roll
C	04144 roll of 2000 labels	yes	70 x 33	40	1 roll
C	04295 roll of 1000 labels	no	70 x 33	76	1 roll
C	04169 roll of 1000 labels	yes	70 x 33	76	1 roll
E	04336 with information on standards and color change, steam indicator left - roll of 1000 labels	yes	70 x 35	40	1 roll
G	04183 steam indicator right - roll of 1000 labels	yes	70 x 35	40	1 roll
H	04182 steam indicator left - roll of 1000 labels	yes	70 x 35	40	1 roll



A




B



C




D




E



F



I



J


SINGLE LABELS - WITH STEAM INDICATOR SUITABLE FOR NON-WOVEN MATERIAL

Blank with steam indicator (pink to brown). Wound outside and duplex self-adhesive.


REF	Description	Tear-off perforation	Format (mm)	Core-Ø (mm)	Qty
B	04250 roll of 1000 labels	yes	70 x 35	40	1 roll
B	04242 roll of 1000 labels	yes	70 x 35	76	1 roll
C	04180 roll of 1000 labels	yes	70 x 35	40	1 roll
D	04273 with information on standards and color change, color bar red - roll of 1000 labels	yes	70 x 35	76	1 roll
D	04275 with information on standards and color change, green color bar - roll of 1000 labels	yes	70 x 35	76	1 roll
D	04276 with information on standards and color change, yellow color bar - roll of 1000 labels	yes	70 x 35	76	1 roll
D	04274 with information on standards and color change, blue color bar - roll of 1000 labels	yes	70 x 35	76	1 roll
E	04280 with information on standards and color change, vapor indicator left - roll of 1000 labels	yes	75 x 35	76	1 roll
F	04272 vapor indicator left - roll of 1000 labels	yes	70 x 35	76	1 roll
F	04328 vapor indicator left - roll of 2000 labels	yes	70 x 35	76	1 roll
H	04311 with information on standards and color change - roll of 1000 labels	yes	100 x 60	76	1 roll
I	04241 roll of 1000 labels	yes	83 x 40	76	1 roll
J	04093 endless labels, wrapped outside, double self-adhesive, also sticks to non-woven materials	no	70 x 80	50	1 roll




A



B




C




D



E



F



G



H



I

SINGLE LABELS - WITH STEAM INDICATOR (PINK TO BROWN)

Wound on the outside and single self-adhesive.

REF	Description	Format (mm)	Core-Ø (mm)	Qty
A	04066 for Webeco, perforated in the middle, with tear-off perforation	89 x 36	40	1 roll
B	04190 suitable for non-woven material - roll of 1000 labels	60 x 31	40	1 roll
C	04094 suitable for non-woven material - roll of 2000 labels	60 x 31	40	1 roll
D	04265 roll of 1000 labels	60 x 33	40	1 roll
E	04266 roll of 2000 labels	60 x 33	40	1 roll
C	04306 for non-woven material, indicator left - roll similar 1000 labels	70 x 33	40	1 roll
F	04198 for non-woven material - roll of 1000 labels	70 x 35	40	1 roll
F	04236 for non-woven material - roll of 2500 labels	70 x 35	40	1 roll
G	04157 for non-woven material - roll of 2500 labels	70 x 35	40	1 roll
H	04256 for non-woven material, indicator left - roll similar 1000 labels	70 x 35	76	1 roll
I	04021 red, for non-woven material - roll of 2000 labels	90 x 35	40	1 roll




SINGLE LABELS - WITH STEAM INDICATOR (PINK TO BROWN)

Wound on the outside and single self-adhesive. With steam indicator.

REF	Description	Tear-off perforation	Format (mm)	Core-Ø (mm)	Qty
A	04299 for non-woven material, with information on standards - roll of 2000 labels	yes	90 x 35	76	1 roll
B	04133L for non-woven material and for longer warmth indicator left - roll of 750 labels	yes	70 x 38	76	1 roll
C	04038 roll of 2000 labels		63 x 40	40	1 roll

ENDLESS LABELS - WITH STEAM INDICATOR (PINK TO BROWN)


Wound outside. Single use.



REF	Description	Self-adhesive	Roll length	Width (mm)	Core-Ø (mm)	Qty
D	04063 yellow	duplex	80	70	50	1 roll
D	04061 green	duplex	80	70	50	1 roll
D	04060 red	duplex	80	70	50	1 roll
E	04085 suitable for non-woven material	duplex	80	70	50	1 roll
E	04145 suitable for non-woven material	duplex	80	70	76	1 roll
F	04181 suitable for non-woven material, with free border on both sides	duplex	80	70	50	1 roll
E	04142 with special adhesive (minimizes soiling of the cutting knife in the label printer)	duplex	80	70	50	1 roll
E	04268 with special adhesive	duplex	80	70	76	1 roll
D	04062 blue bar	duplex	80	70	50	1 roll
G	04098 continuous label on reel with steam indicator	single	100	60	76	1 roll




A




B




C



D




E



F



G



H

LABELS - WITH FO/EO AND STEAM INDICATOR

REF	Description	Format (mm)	Core-Ø (mm)	Qty
A	04185 single labels, wound inside, double self-adhesive - roll of 1000 labels	60 x 31	40	1 roll
B	04065 individual labels for Webeco, wound outside, perforated in the middle, with tear-off perforation, with back slit - roll of 1000 labels	86 x 36	40	1 roll
C	04045 continuous labels (EO = white, DI = pink), wound outside, double self-adhesive	70	50	1 roll
D	04046 continuous labels (FO= violet, DI = pink), wound outside, double self-adhesive	70	50	1 roll

LABELS WITH H₂O₂ AND STEAM INDICATOR


With tear-off perforation, for single use, folded self-adhesive.

REF	Description	Format (mm)	Core-Ø (mm)	Qty
E	04135 singel labels, wound outside	70 x 35	40	1 roll
F	04279 singel Tyvek labels, wound inside	70 x 35	76	1 roll
G	04097 continuous label, wound outside	70	50	1 roll


OTHER SINGLE LABELS - WITH PUNCHED EDGE

With steam indicator (pink to brown), zigzag fold, with tear-off perforation. Single use.

REF	Description	Format (mm)	Qty
H	04004 duplex self-adhesive, universal adhesive, printing: content, name, sterile date and patient.	99 x 37	2500 pcs



A



A



C

Foil labels allow you to identify your containers and sieve tray. They offer durability and adhesion to smooth surfaces, even after the cleaning and disinfection process. Key Surgical offers you different versions in numerous sizes compatible with the most common label printers.

FOIL LABELS - BLANK - WITHOUT INDICATOR

Heat-resistant up to 134 °C, suitable for washer-disinfectors, wound outside, single self-adhesive.


REF	Description	Tear-off perforation	Format (mm)	Core-Ø (mm)	Qty
A 04255	yellow, suitable for washer-disinfectors and dryers, strong adhesive	yes	60 x 25	40	1 roll
A 04222	green, suitable for washer-disinfectors, strong adhesive	yes	60 x 25	40	1 roll
A 04221	red, not washable, strong adhesive	yes	60 x 25	40	1 roll
B 04271	white, suitable for washer-disinfectors, strong adhesive	yes	50 x 23	25	1 roll
B 04176	white, suitable for washer-disinfectors, strong adhesive	yes	60 x 25	40	1 roll
B 04218	suitable for washer-disinfectors	no	40 x 28	25	1 roll
B 04174	suitable for washer-disinfectors and dryers, strong adhesive	yes	40 x 28	40	1 roll
B 04186	suitable for washer-disinfectors, with one-sided barrel flap, can be removed without leaving any residue		55 x 35	50	1 roll
C 04227	white, suitable for washer-disinfectors, strong adhesive	yes	50 x 23	40	1 roll

FOIL LABELS FOR IDENTIFYING AND SEALING CONTAINERS

Can be removed without any residue, single self-adhesive, wound outside. Single use.

REF	Description	Format (mm)	Core-Ø (mm)	Qty
01001	with steam indicator	102 x 25	76	1 roll
01002	without steam indicator	102 x 25	76	1 roll
04096	without steam indicator	75 x 35	76	1 roll

A



A THERMO TRANSFER FOILS


Thermal transfer reels with various coatings and features, suitable for your printer and labels. We recommend the premium resin which is particularly suitable for film carrier material (except 09858, which is only for paper). Wax/resin features are suitable for the most common requirements of the central sterile services department (CSSD) and especially for paper labels, and the wax feature is suitable for simple applications without thermal requirements. Coated on the outside or inside, so you are guaranteed to find the right foil reel for your printer in our range.

REF	Coating	Quality	Format (mm)	Sleeve length (mm)	Sleeve-Ø (mm)	Autoclave resistant	Qty
09476	outside	resin	85 x 74	110	12.5	up to 150 °C	1 roll
09213	outside	wax	85 x 74	110	12.5	yes	1 roll
09452	outside	wax/resin	85 x 74	110	12.5	no	1 roll
09212	outside	wax	110 x 74	110	25	yes	1 roll
09853	outside	wax/resin	110 x 74	110	12.5	no	1 roll
09080	outside	wax	70 x 300	70	25	yes	1 roll
09270	outside	wax/resin	70 x 300	70	25	yes	1 roll
09069	inside	wax	76 x 300	76	25	yes	1 roll
09852	inside	wax/resin	76 x 300	76	25	no	1 roll
09486	outside	resin	80 x 300	80	25	up to 150 °C	1 roll
09077	outside	wax	80 x 300	80	25	yes	1 roll
09854	outside	wax/resin	80 x 300	80	25	no	1 roll

ADHESIVE LABELS

Adhesive labels come in a DIN format. They are single self-adhesive.

REF	Description	Indicator	Format (mm)	Qty
A 04016	for laser printers, single self-adhesive - 25 x 44-label DIN A4 sheets	steam	48 x 25	1100 pcs
A 04017	for laser printers, single self-adhesive - 25 x 21-label DIN A4 sheets	steam	64 x 38	525 pcs



A

PRE-PRINTED LABELS

Pre-printed labels provide visual support for communication between the individual departments. Additionally, warning signs ensure the protection of your staff, for example for burns on hot surfaces, by changing their color when they get hot. Labels for disinfectant bottles help you to comply with labelling requirements.

REF	Description	Imprint	Color	Format (mm)	Roll width (mm)	Core-Ø (mm)	Qty	
A	04326	wound outside, made of PE film	Anbruch am:____ Verwendbar bis:____	yellow	70 x 30	102	40	1 roll
B	04281	wound inside, single self-adhesive serially numbered - roll of 500 labels	Sicherheitssiegel / Security Seal	green	90 x 30	102	76	1 roll
C	04240	wound inside, single self-adhesive serially numbered - roll of 500 labels	desinfiziert, wenn Etikett intakt	red	90 x 30	102	76	1 roll
D	04156	wound outside, single self-adhesive- roll of 1000 labels	gereinigt + desinfiziert	blue	100 x 40	102	40	1 roll
E	04238	wound outside, single self-adhesive, heat-resistant up to 134 °C - roll of 1000 labels	unvollständig	orange	100 x 40	102	40	1 roll
F	04288	wound outside, single self-adhesive roll of 1000 labels	kontaminiert	red	100 x 40	102	40	1 roll





A



C




D



E

PRE-PRINTED LABELS

REF	Description	Imprint	Color	Format (mm)	Roll width (mm)	Core-Ø (mm)	Qty
A 04215	wound inside, single self-adhesive, sterilizable up to 134 °C - roll of 1000 labels	loan sieve	yellow	100 x 40	104	40	1 roll
C 04228	wound inside, single self-adhesive, not sterilizable - roll of 1000 labels	defect	red	40 x 50	43	76	1 roll
04229	wound inside, single self-adhesive, not sterilizable - roll of 1000 labels	not sterile	red	40 x 50	43	76	1 roll
D 04177	wound inside, single self-adhesive, not sterilizable - roll of 1000 labels	thermal disinfected	red	40 x 50	43	40	1 roll
E 04332	self-adhesive - removable adhesive, 2 lines	sterile mp	green	150 x 90	NA		250 pcs




A




B



C



D



E

SINGLE LABELS FOR LABELLING DEVICES

REF	Description	Color	Format (mm)	Qty
A 04034	for labelling device REF 04011, with pre-printed text: 'Sterilized on/usable until'	blue	37 x 19	1 roll
B 04072	METO - for labelling device 04019, with steam indicator	blue	29 x 28	12 rolls
B 04073	METO - for labelling device 04019, with steam indicator	yellow	29 x 28	12 rolls
B 04071	METO - for labelling device 04019, with steam indicator	green	29 x 28	12 rolls
B 04068	METO - for labelling device 04019, with steam indicator	orange	29 x 28	12 rolls
B 04070	METO - for labelling device 04019, with steam indicator	red	29 x 28	12 rolls
B 04069	METO - for labelling device 04019, with steam indicator	white	29 x 28	12 rolls
C 04054	for labelling device REF 04014, without indicator	blank	23 x 31	12 rolls
D 04041	for labelling device REF 04014, without indicator	blue	23 x 31	12 rolls
E 04052	for labelling device REF 04014, with steam indicator	blue	23 x 31	12 rolls
E 04053	for labelling device REF 04014, with steam indicator	yellow	23 x 31	12 rolls
D 04040	for labelling device REF 04014, without indicator	green	23 x 31	12 rolls
E 04051	for labelling device REF 04014, with steam indicator	green	23 x 31	12 rolls
E 04048	for labelling device REF 04014, with steam indicator	orange	23 x 31	12 rolls
E 04050	for labelling device REF 04014, with steam indicator	red	23 x 31	12 rolls
E 04049	for labelling device REF 04014, without indicator	white	23 x 31	12 rolls

LABELLERS AND ACCESSORIES

REF	Description	Qty
A	04018 MONARCH - ink refill for "Monarch" labeller REF 04014	2 pcs
B	04232 ink refill for labeller REF 04230 – 3 lines	1 pc
C	04119 METO ink refill for "Meto" labeller REF 04019	2 pc
D	04014 MONARCH "Monarch" labeller – 3 lines	1 pc
E	04019 METO "Meto" labeller – 3 lines (12 characters per line (figures only); optionally available with text and special characters)	1 pc
F	04231 labeller for manual marking – 2 lines (10 characters per line (figures only)	1 pc
G	04230 labeller for manual marking – 3 lines (11 characters per line (figures only)	1 pc



A



B



C



D




E



F



G





MARKERS

REF	Description	Ink color	Line thickness (mm)	Qty
A 05504	edding 8400 permanent marker, lightfast	black	0.5 – 1	10 pcs
05513	permanent marker, sterilization-proof	red	0.75	10 pcs
B 05505	permanent marker, sterilization-proof	black	0.75	10 pcs
C 09089	edding 8010 X-ray marker	white	0.8	10 pcs
D 09087	edding 400 permanent marker, lightfast	blue	1	10 pcs
E 09170	edding 400 permanent marker, lightfast	red	1	10 pcs
F 09086	edding 400 permanent marker, lightfast	black	1	10 pcs
G 09091	edding 8020 skin marker, with round tip, not water resistant	blue	1	10 pcs
H 09090	edding 8020 skin marker, with round tip, not water resistant	red	1	10 pcs
I 05514	SteriLiner permanent marker, sterilization-proof	black	1	10 pcs
J 05515	SteriLiner XL permanent marker, with chisel tip, sterilization-proof	black	2 – 5	10 pcs
K RT-209P	Sharpie lab marker	black		12 pcs

BROTHER P-TOUCH

REF	Model	Description	Qty
M 09500	PT-D600VP	<ul style="list-style-type: none"> • German QWERTZ keyboard • color display, 3 lines x 20 characters • up to 24mm tape width • 280 displayable characters, phrase memory with 2800 characters • automatic saving of the last entry • 14 fonts, 10 print styles, 6 - 48 point size, 99 frame types and background patterns • vertical print, multiple printouts up to 99 times, multiline printouts up to 7 lines, mirror writing, tabulator function • 25 pre-programmed label templates (+ downloadable template), adjustable label length, 9 standard barcodes, multi-block printing up to 5 blocks • automatic cutter • recognizes the color of the inserted tape • USB Interface Layout design on the PC • logo and symbol printing in up to 180 x 360 dpi 	1 pc



LABELLING TAPES FOR THE P-TOUCH LABELLING MACHINE

Heat-resistant up to 134 °C, washable, only suitable for steam sterilization. Run length 8 meters.

REF	Tape width (mm)	Laminated	Ribbon color	Font color	Qty
A 09511	18	yes	white	black	1 pc
B 09514	18	yes	yellow	black	1 pc
C 09544	18	no	white	black	1 pc
D 09543	24	no	white	black	1 pc
E 09502	24	yes	white	black	1 pc
F 09503	24	yes	red	black	1 pc
G 09505	24	yes	yellow	black	1 pc
H 09506	24	yes	green	black	1 pc


STENCIL TAPES FOR FOR CHEMICAL ETCHING

REF	Font color	Ribbon color	Format (mm)	Run length (m)	Qty
J 09536	black	white	24	3	1 pc

STATUS DOCUMENTATION

With tear-off perforation, for single use, double self-adhesive.

REF	Description	Format (mm)	Qty
K 09072	tray documentation block with steam indicator - block of 100 sheets	100 x 100	1 block
L 09071	tray documentation pad without steam indicator made from special medical paper - pad of 100 sheets	100 x 100	1 block
M 09214	tray control card with steam indicator, reverse side printing: TIME OUT	120 x 75	1000 pcs



J

K

L

M

