

## INSTRUMENT PROTECTION

The protection of instruments has never been more important. As a result of increasing centralization, the merger of CSSDs and external reprocessing, transport routes are becoming ever longer. This means that sensitive and complex instruments have to be protected not only within your own hospital, but often when being transported off-site.

Instrument Sleeves are made of special cardboard or transparent film and are used to protect high-quality instruments from damage during preparation, transport and storage. The shape of the protectors is not directly adapted to the respective instrument. This has the advantage that these protective sleeves are suitable for the most common instruments as well as for special shapes (such as hooks or knives).

Instrument Tip Caps are designed to extend the life of your instruments. Our tip caps protect your valuable instruments from unnecessary wear or damage, helping reduce repair costs while giving you the peace of mind that your instrument investment is being protected. In addition to protecting and extending the life of your valuable instruments, tip caps also help to protect sterilization packs from puncture and protect employees from sharp instruments.

Available in a variety of shapes and sizes to fit all your instrument protection needs, we have the perfect tip cap for each of your various instruments, from flat caps to twin caps to retractor caps, and they are available with or without perforations, depending on your preference. Our tip caps are radiopaque, single use, and compatible with steam sterilization.

[Instrument Sleeves](#)

[Instrument Protection Caps Round](#)

[Instrument Organizers](#)

[Instrument Protection Caps Flat](#)

[Special Instrument Protection Caps](#)

Unless otherwise stated, our products are not made with natural rubber latex.



A



B



C



D




## Instrument Protection – Instrument Sleeves

### INSTRUMENT SLEEVES

Instrument Sleeves are used to protect high-quality instruments from damage during preparation, transport, and storage. They are therefore suitable for the most common instruments as well as for special shapes (such as hooks or knives). For single use.

REF	Description	Material	Format* (mm)	Qty
A	01455 hybrid instrument sleeve	cardboard/ film	120 x 32	100 pcs
B	01457 hybrid instrument sleeve to hold hinged instruments open	cardboard/ film	160 x 40	100 pcs
C	01476 instrument sleeve	film	60 x 10	1000 pcs
D	01471 folding borders can be folded inwards to prevent narrow instruments from slipping out	cardboard	105 x 25	1000 pcs
	01477 folding borders can be folded inwards to prevent narrow instruments from slipping out	film	140 x 50	100 pcs
E	01479 for amputation knives, with 8 holes	cardboard	123 x 33	100 pcs
F	01478 for amputation knives, with 8 holes	cardboard	246 x 33	100 pcs
	PPT001 white card		160 x 50	2000 pcs
	PPT002 white card		160 x 50	1000 pcs
	PPT004 white card		160 x 50	1000 pcs

\* Useful length x max. instrument width



K

Instrument Protection

## Instrument Protection – Round Caps



Guards and caps are used to protect premium instruments from damage during transportation, handling and sterilization. Furthermore, they help to protect personnel from cut wounds and minimize reprocessing costs caused by rips or punctures in sterilization packaging.

### ROUND CAPS

For single use.

REF	Description	Color	Format (mm) Ø	Qty
A	3-04-01 perforated	clear	1.6 x 19	100 pcs
B	3-08-01 perforated	white transparent	1.6 x 19	100 pcs
C	3-15-01 perforated	white	1.6 x 19	100 pcs
D	01856 for Verres needles	clear	2.0 x 19	50 pcs
E	3-04-02 perforated	clear	2.0 x 19	100 pcs
F	3-08-02 perforated	blue transparent	2.0 x 19	100 pcs
G	3-15-02 perforated	blue	2.0 x 19	100 pcs
H	3-04-03 perforated	clear	2.8 x 19	100 pcs
I	3-08-03 perforated	green transparent	2.8 x 19	100 pcs
J	3-15-03 perforated	green	2.8 x 19	100 pcs
K	3-04-04 perforated	clear	3.2 x 25	100 pcs
L	3-08-04 perforated	red transparent	3.2 x 25	100 pcs
M	3-15-04 perforated	red	3.2 x 25	100 pcs
N	3-04-99 perforated	clear	4.7 x 38	100 pcs
O	3-04-05 perforated	clear	5.0 x 25	100 pcs
P	3-08-05 perforated	yellow transparent	5.0 x 25	100 pcs
Q	3-15-05 perforated	yellow	5.0 x 25	100 pcs
R	3-04-06 perforated	clear	10 x 19	100 pcs
S	3-08-06 perforated	purple transparent	10 x 19	100 pcs
T	3-15-06 perforated	purple	10 x 19	100 pcs
U	01859 for obturators	clear	10 x 20	25 pcs
V	3-07-03 for MIS instruments	green	4.3 x 29	50 pcs
W	3-07-05 for MIS instruments	yellow	4.3 x 51	50 pcs
X	01857 for obturators	clear	5.0 x 25	25 pcs
Y	01858 for obturators	clear	5.0 x 43	20 pcs
Z	01860 for obturators	clear	9.1 x 66	20 pcs
AA	3-07-09 for MIS instruments	clear	9.0 x 66	50 pcs

PRODUCTS SHOWN ACTUAL SIZE

## Instrument Protection – Tapered Caps

### ROUND CAPS – TAPERED CAPS

For single use.

REF	Color	Format (mm) Ø	Qty
A	white transparent	2.0 - 3.0 x 38	70 pcs
B	white transparent	2.0 - 3.0 x 51	70 pcs
C	blue transparent	3.0 - 4.0 x 38	70 pcs
D	blue transparent	3.0 - 4.0 x 51	70 pcs
E	green transparent	4.0 - 5.0 x 51	70 pcs
F	green transparent	4.0 - 5.0 x 38	70 pcs
G	red transparent	5.0 - 6.0 x 51	70 pcs
H	red transparent	5.0 - 6.0 x 38	70 pcs
I	yellow transparent	6.0 - 7.0 x 38	70 pcs
J	yellow transparent	6.0 - 7.0 x 51	70 pcs
K	purple transparent	7.0 - 8.0 x 38	70 pcs
L	purple transparent	7.0 - 8.0 x 51	70 pcs
M	orange transparent	8.0 - 9.0 x 38	70 pcs
N	orange transparent	8.0 - 9.0 x 51	70 pcs
O	brown transparent	9.0 - 10 x 38	70 pcs
P	brown transparent	9.0 - 10 x 51	70 pcs
Q	black transparent	10 - 11 x 38	70 pcs
R	black transparent	10 - 11 x 51	70 pcs
S	green transparent	13 - 14 x 38	50 pcs
T	green transparent	15 - 16 x 51	50 pcs
U	yellow transparent	15 - 16 x 38	50 pcs
V	yellow transparent	19 - 20 x 51	50 pcs
W	orange transparent	15 - 16 x 38	50 pcs
X	orange transparent	19 - 20 x 51	50 pcs



PRODUCTS SHOWN ACTUAL SIZE

## Instrument Protection – Flat Caps



### FLAT CAPS

For single use.

REF	Description	Color	Format (mm)	Qty
A	3-04-07 for osteotomes	clear	2.0 x 9.0 x 25	100 pcs
B	3-08-07 for osteotomes	orange transparent	2.0 x 9.0 x 25	100 pcs
C	3-15-07 for osteotomes	orange	2.0 x 9.0 x 25	100 pcs
D	3-04-08 for osteotomes	clear	2.0 x 16 x 25	100 pcs
E	3-08-08 for osteotomes	brown transparent	2.0 x 16 x 25	100 pcs
F	3-15-08 for osteotomes	brown	2.0 x 16 x 25	100 pcs
G	3-04-09 for osteotomes	clear	3.0 x 25 x 25	100 pcs
H	3-08-09 for osteotomes	black transparent	3.0 x 25 x 25	100 pcs
I	3-15-09 for osteotomes	black	3.0 x 25 x 25	100 pcs
J	3-04-11 for osteotomes	clear	6.4 x 19	50 pcs
K	3-08-11 for osteotomes	white transparent	6.4 x 19	100 pcs
L	01807 for osteotomes	clear	9.5 x 25	50 pcs
M	3-08-12 for osteotomes	blue transparent	9.5 x 25	100 pcs
N	01808 for osteotomes	clear	13 x 25	50 pcs
O	3-08-13 for osteotomes	green	13 x 25	100 pcs
P	3-08-15 for osteotomes	yellow	18 x 25	100 pcs
Q	3-04-15 for osteotomes	clear	19 x 25	50 pcs
R	3-04-09 for osteotomes	clear	25 x 25	50 pcs
S	3-04-17 for osteotomes	clear	32 x 32	50 pcs
T	3-08-17 for osteotomes	purple transparent	32 x 32	100 pcs
U	3-04-18 for osteotomes	clear	38 x 38	50 pcs
V	3-08-18 for osteotomes	red	38 x 38	100 pcs
W	3-04-19 for osteotomes	clear	45 x 32	100 pcs
X	01813 for osteotomes	clear	45 x 32	50 pcs
Y	3-08-19 for osteotomes	orange	45 x 32	100 pcs

PRODUCTS SHOWN ACTUAL SIZE

## Instrument Protection – Skin Hook Caps and Retractor Caps



PRODUCTS SHOWN ACTUAL SIZE

### SKIN HOOK CAPS

Protect delicate single and double prong skin hook instrumentation. Radiopaque. Compatible with sterilization. Single use.

REF	Description	Color	Format (mm)	Qty
A	3-10-03 double prong. 2mm apart	clear	4.0 x 19	100 pcs
B	3-10-01 single prong. 2.5-3.5mm hook	clear	4.8 x 19	100 pcs
C	3-10-02 single prong. 5mm hook	clear	6.3 x 19	100 pcs
D	3-10-04 double prong. 5mm apart	clear	6.7 x 19	100 pcs
E	3-10-05 double prong. 7mm apart	clear	9.1 x 19	100 pcs

### SELF-RETAINING AND GELPI RETRACTOR CAPS


Protect a variety of self-retaining retractors including Weitlaner, D'Errico-Adson, Beckman, and Adson retractors. Radiopaque. Compatible with steam sterilization. Single use.

REF	Description	Color	Format (mm)	Qty
F	3-09-07 self-retaining retractor cap	blue transparent	22 x 20	100 pcs
G	3-09-08 self-retaining retractor cap	purple	25 x 25	100 pcs
H	3-09-09 self-retaining retractor cap	orange	25 x 32	100 pcs
I	3-09-10 self-retaining retractor cap	green	45 x 13	100 pcs
J	3-09-03 Gelpi retractors	clear	51	100 pcs
K	3-09-04 Gelpi retractors	clear	76	100 pcs
L	3-09-05 Gelpi retractors	clear	38	100 pcs

## Instrument Protection – Twin Caps

K

Instrument Protection



### TWIN CAPS

For single use.

REF	Description	Color	Format (mm) Ø	Qty
A	3-06-01 NA	white transparent	1.6 x 19	50 pcs
B	3-06-02 NA	blue transparent	2.0 x 19	50 pcs
C	3-06-03 NA	green transparent	2.8 x 19	50 pcs
D	3-06-07 NA	orange transparent	2.0 x 25	50 pcs
E	3-06-04 NA	red transparent	3.2 x 25	50 pcs
F	3-06-05 NA	yellow transparent	5.0 x 25	50 pcs
3-06-50	assortment (1 package of each) 3-06-01, 3-06-02, 3-06-03, 3-06-04, 3-06-05, 3-06-07	clear	-	1 set


### AIRFLOW CAPS™

An alternative to perforated caps, grips on the inside of the cap allow more airflow during sterilization process and assist in easy removal. Radiopaque. Compatible with steam sterilization. Single use.

REF	Color	Format (mm)	Qty
G	clear	2.0 x 16 x 25	100 pcs
H	clear	4.8 x 25	100 pcs
I	clear	9.5 x 25	100 pcs
J	clear	10 x 20	100 pcs
K	clear	25 x 25	100 pcs

PRODUCTS SHOWN ACTUAL SIZE

## Instrument Protection – Instrument Organizers



### INSTRUMENT ORGANIZERS

Designed to help with efficient organization and protection of surgical instruments during the sterilization process. Must be used in conjunction with a sterile barrier system, the organizers feature multiple slots and rows to hold and protect your instruments.

REF	Description	Length (mm)	Width (mm)	Qty
<b>A</b>	CS-DH100 2 rows of 6 pockets	317.5	438.1	100 pcs
<b>B</b>	CS-DH102 2 rows of 6 pockets with slots in pocket	304.8	444.5	100 pcs
<b>C</b>	CS-DH200 2 rows of 6 notched pockets with tape closure strips	273	482.6	100 pcs
<b>D</b>	CS-DH300 2 rows of 6 pockets with foldover flap	412.7	444.5	100 pcs
<b>E</b>	CS-DH505 single pocket	127	146	100 pcs
<b>F</b>	CS-DH512 single pocket	127	305	100 pcs

