

PERSONAL PROTECTIVE EQUIPMENT

Face Shields are available in several styles to fit the preferences of every individual. We offer full and short face shields, with masks, with Tyvek™ expanders, with ties and with ear loops. All of our face shields offer different types of protection. Barrier Gowns are lightweight, comfortable, and breathable.

Nitrile gloves are powder-free, available in 24cm, 30cm and 40cm. They are treated with an antibacterial and anti-mildew agent.

We also offer sterilizer mitts made in silicone. These provide you with a protection barrier for removing hot trays and racks.

Protective Apparel

Face Shields and Masks

Gloves



Unless otherwise stated, our products are not made with natural rubber latex.



Personal Protective Equipment – Protective Apparel

CONVENIENT GOWN AND PROTECTIVE / VISITOR GOWN

We have a wide range of disposable gowns in various materials designed to protect personnel against contact with cleaning materials of weak action or prolonged contact with water. Our range is tailored to your needs, whether coated fronts and sleeves, for wet areas, convenient disposable visitor gowns or for short and less critical applications. We offer a large selection in all common sizes, as well as special XXL gowns. Single use.

REF	Description	Material	Color	Size (cm)	Weight (gms)	Qty
A 09437	spunbonded gown, non sterile, for single use, breathable, water resistant, long sleeves with elastic cuffs	mixed fleece	blue	130 long	72	50 pcs
EN10120B	uncoated, with thumb loops	SMS	blue	198 x 105	30	100 pcs
B EN10120G	coated, with thumb loops	SMS	yellow	198 x 105	30	100 pcs
B EN10120GL	coated, with thumb loops	SMS	yellow	202 x 115	30	50 pcs
EN10211.0	with knitted cuffs, completely coated gown	fleece	yellow	136 x 140	31	100 pcs
EN10120GB	coated, with knitted cuffs	SMS	yellow	198 x 105	30	100 pcs

CONVENIENT GOWN

For special processes, we offer different examination gowns. The sleeves are made of woven material, which offers more protection, while the gown's back is made of a breathable non-woven material to provide comfort. Non-sterile. Single use.

REF	Description	Color	Size (cm)	Weight (gms)	Qty
EN10120.1G	with elastic cuffs	yellow	136 x 140	26	100 pcs
EN10120.0PLUS	with knitted cuffs	yellow	136 x 140	36	100 pcs
EN10120.3	with knitted cuffs	yellow	140 x 150	26	100 pcs
EN10120.0XL	with knitted cuffs	yellow	160 x 160	26	100 pcs





Personal Protective Equipment – Protective Apparel

EXAMINATION GOWNS

Our examination gowns are equipped with cuffs and a color coded tie belt. They also have inside reinforcements at the seams in the chest and arm areas to reduce the risk of contamination. They are made from SMS and have a material thickness of 35 g/m². Single use.

REF	Color	Size (cm)	Weight (gms)	Qty
A EN10119M	light blue	120 x 140	35	50 pcs
A EN10119L	light blue	125 x 155	35	50 pcs
A EN10119XL	light blue	130 x 160	35	50 pcs

PROTECTIVE / VISITOR GOWNS

Our gowns are ideal for hospital visitors. They are made of fleece with elastic or knitted cuffs. Single use.

REF	Description	Color	Size (cm)	Weight (gms)	Qty
B EN10500.0	with knitted cuffs	green	137 x 140	20	100 pcs
EN10500.1	with elastic cuffs	green	136 x 140	20	100 pcs
EN10500.2	with knitted cuffs	light blue	140 x 138	20	100 pcs
B EN10525.0	with knitted cuffs	green	136 x 140	25	100 pcs
EN10525.1	with elastic cuffs	green	136 x 140	25	100 pcs
B EN10530.0	with knitted cuffs	green	136 x 140	30	100 pcs
EN10530.1	with elastic cuffs	green	140 x 140	30	100 pcs
EN10535.1	with elastic cuffs	green	136 x 140	35	100 pcs

APRONS AND GASTROSCOPY BIBS

REF	Description	Size (cm)	Qty
EN10123B	colonoscopy aprons, made of PE film	81 x 125	100 pcs
EN10124	castroscopy bib, made of PE film	37 x 66	500 pcs

Personal Protective Equipment – Face Shields & Masks



FACE PROTECTION AND MASKS

Face and mouth protection are essential, especially when working in reprocessing areas. We are committed to high quality products such as full face visors, protective goggles, and face masks. Our face shields are designed to protect personnel against contact with cleaning materials or prolonged contact with water.

MASKS

Our FFP2 mouth and nose protection masks offer additional safety on a daily basis.

REF	Description	Qty
A 09441	mask with visor, blue-white with elastic ear-loop	25 pcs
B 09440	mask, blue-white, with elastic ear-loop	50 pcs



FULL FACE MASKS

REF	Description	Material	Qty
C 05TD3120ST	clear vision comfort face shield eco long, two face protective layer, gray foam, 21.5cm long	PTE	100 pcs
C EN10152AF	full face visor with foam pads and elastic strap	PTE	200 pcs

FACE VISORS

REF	Description	Color	Qty
D 09081	frame for 2-piece safety goggles	yellow	10 pcs
E 09431	frame for 2-piece safety goggles	pink	10 pcs
F 09082	visors for 2-piece safety goggles made of plastic	NA	25 pcs
05TD0204	frame for 2-piece safety goggles suitable for REF 05TD0161	multicolor	50 pcs



Personal Protective Equipment – Face Shields & Masks

FACE PROTECTION SET- CLEARVIZOR®

Set consists of:

- 1x frame
- 10x foil shields
- No distortion of the field of vision
- Offers health protection during work
- Anti-glare
- Anti-static
- Equipped with anti-fog properties
- Scratch-resistant
- Resistant to liquids and chemicals
- Safe (resistant to breaking and chipping)
- The visor is protected on both sides with adhesive foils (incl. removal aids)
- Light and ergonomic

REF

- A** 09714
- B** 09715
- C** 09716
- D** 09717
- E** 09718
- F** 09719

Color

- white
- blue
- transparent
- pink
- yellow
- green

Qty

- 1 set
- 1 set
- 1 set
- 1 set
- 1 set
- 1 set



A



Personal Protective Equipment – Gloves

PROTECTIVE GLOVES

Our protective gloves range satisfies the requirements of various working areas. Our silicone gloves meet high quality standards and have been optimized for the tough demands of the CSSD in terms of ergonomics and performance. Flexibility, resistance, and comfort are our key quality criteria for heat-protective gloves.

HEAT PROTECTION GLOVES

REF	Size	Length (cm)	Material	Heat resistance	W/D compatible	Color	Qty
A 09368	universal	35	silicone	-40 °C up to +240 °C	Yes	light blue	1 pair
B 09263	universal	37	silicone	-60 °C up to +260 °C	Yes	blue	1 pair

B



K

Personal Protective Equipment

Personal Protective Equipment – Gloves

A



A 24CM EXAMINATION GLOVES

Our examination gloves are made of nitrile, have a standard length of 24cm, and are available in different sizes. Single use.

REF	Size	Color	Qty
RNB01S	S (6½ – 7)	blue	200 pcs
RNB01M	M (7½ – 8)	blue	200 pcs
RNB01L	L (8½ – 9)	blue	200 pcs
RNB01XL	XL (9½ – 10)	blue	180 pcs

B 30CM EXAMINATION GLOVES

Our protective gloves are made of nitrile, have a standard length of 30cm, and are available in different sizes. Single use.

REF	Size	Color	Qty
MED1100.90	S (6½ – 7)	blue	100 pcs
MED1100.93	XL (9½ – 10)	blue	100 pcs

40CM EXAMINATION GLOVES

Our protective gloves are made of nitrile, have a standard length of 40cm, and are available in different sizes. Single use.

REF	Size	Color	Qty
09472	S (6½ – 7)	blue	50 pcs
09473	M (7½ – 8)	blue	50 pcs
09474	L (8½ – 9)	blue	50 pcs
09475	XL (9½ – 10)	blue	50 pcs

CUFFS

REF	Description	Material	Size	Color	Qty
EN10153	disposable cuffs	PE film	universal	blue	100 pcs

B



K

Personal Protective Equipment