

WRAPS & LINERS

Our flexible production methods ensure that the best sizes of wrap can be produced for our Customers, even if these are not industry standard sizes, without compromising performance or protection levels.

Our Barrier Wraps provide a bacterial barrier and protect the contents of a set or pack from contamination. When classed as alcohol repellent, barrier wraps can be used within the sterile field.

Choice Wrap is designed to reduce the risk of rips, tears, and wet packs. Choice SMS are made of polypropylene. Choice Wraps can be interleaved in blue or green to meet the sterilization service user requirements.

All Transport Wraps must be used in conjunction with one compliant Barrier Wrap to completely prevent the ingress of microorganisms. Use of a single barrier layer with Transport Wrap can be a cost-effective solution that reduces damage in handling and transport. If adopted, changes to both wrapping and unwrapping procedures may be required to ensure aseptic opening.

Choice Wrap

Double Wrap


Crepe Wrap

Tray Liners

Indicator Tape and Dispenser

Transport Wrap

Unless otherwise stated, our products are not made with natural rubber latex.





CHOICE WRAP 55GSM

A fully synthetic SMS product designed to offer flexibility and enable users to pack conventionally or with 'two layers' at once. Choice Wrap is strong and available with alcohol repellence in both green and blue for use within the sterile field. Wraps are also available with blue and green interleaved. Compatible with steam sterilization. Single use.

A CHOICE WRAP BLUE

REF	Format (cm)	Qty
CB060	60 x 60	350 pcs
CB075	75 x 75	300 pcs
CB090	90 x 90	300 pcs
CB100	100 x 100	200 pcs
CB120	120 x 120	120 pcs
CB140	140 x 140	100 pcs
CB1315	130 x 150	100 pcs

B CHOICE WRAP GREEN

REF	Format (cm)	Qty
CG060	60 x 60	350 pcs
CG075	75 x 75	300 pcs
CG090	90 x 90	300 pcs
CG100	100 x 100	200 pcs
CG120	120 x 120	120 pcs

CHOICE WRAP COMBINATION (GREEN /BLUE)


REF	Format (cm)	Qty
CG/CB060	60 x 60	500 pcs
CG/CB075	75 x 75	200 pcs
CG/CB090	90 x 90	200 pcs
CG/CB100	100 x 100	200 pcs
CG/CB120	120 x 120	100 pcs

A



B





Wraps & Liners – Choice Wrap



A CHOICE BLUE ULTRA 73GSM

Choice Blue Ultra is extra strong, drapeable, and available with alcohol repellence for use within the sterile field. Compatible with steam sterilization. Single use.

REF	Format (cm)	Qty
CBU120	120 X 120	120 pcs
CBU140	140 x 140	100 pcs
CBU150	150 x 150	100 pcs
CBU1014	100 x 140	120 pcs

B CHOICE GREEN WRAP 55GSM


A fully synthetic SMS product designed to offer flexibility and enable users to pack conventionally or with 'two layers' at once. Available with alcohol non repellent in green for use within the sterile field. Compatible with steam and EO sterilization. Single use.

REF	Format (cm)	Qty
CGX090	90 x 90cm	300 pcs
CGX0100	100 x 100cm	200 pcs

A



B





Wraps & Liners – Double Wrap Choice Together



A RELIANCE® SOLO WRAP

A two-sheet wrap system, with two colors joined together along the sides. Each sheet is a five-layer nonwoven tray wrap, manufactured from spun-bonded and melt blown polypropylene, providing a wet and dry microbiological barrier. Bonding the wrap simplifies the wrapping process whilst increasing durability to avoid rips, tears or punctures. Wraps are also available interleaved with other compliant barrier layers, or transportation wrap. Compatible with steam and EO sterilization. Single use. Choice Wrap variants meet ISO 11607-1 and EN 868-2 requirements.

RELIANCE® SOLO WRAP 200B/G

REF	Format (cm)	GSM	Qty
CT43LGT-90	90 x 90	43gsm	50 pcs
CT43LGT-100	100 x 100	43gsm	50 pcs
CT43LGT-140	140 x 140	43gsm	25 pcs


RELIANCE® SOLO WRAP 300B/G

REF	Format (cm)	GSM	Qty
CT50MED-75U	75 x 75	50gsm	66 pcs
CT50MED-90U	90 x 90	50gsm	54 pcs
CT50MED-120U	120 x 120	50gsm	54 pcs
CT50MED-140U	140 x 140	50gsm	36 pcs

RELIANCE® SOLO WRAP 400B/G

REF	Format (cm)	GSM	Qty
CT60HVY-120U	120 x 120	60gsm	36 pcs
CT60HVY-140U	140 x 140	60gsm	36 pcs

A





RELIANCE® NON-WOVEN WRAP 60GSM

Non-woven wrap is high-performance sterilization material that is suitable for inner and outer wrapping of medium to large trays and packs. Wet laid wrap with 50% cellulose and 50% synthetic fibers that improve strength and increase softness. These products are also available in interleaved green and blue wrap. Single use.

REF	Color	Format (cm)	Qty
A RELBLU100U	blue	100 x 100	204 pcs
RELBLU120	blue	120 x 120	100 pcs
B RELGRN100U	green	100 x 100	204 pcs
RELGRN1315U	green	130 x 150	120 pcs

RELIANCE® CREPE WRAP 60GSM

Crepe wrap is suitable for lighter instruments and sets. Manufactured from 100% medical quality cellulose fibers and combines most economical wrapping with good drapeability. Compatible with steam and EO sterilization. Single use.

REF	Color	Format (cm)	Qty
C GRNCREPE060N	green	60 x 60	750 pcs
GRNCREPE075NU	green	75 x 75	174 pcs
GRNCREPE090U	green	90 x 90	126 pcs
GRNCREPE100U	green	100 x 100	126 pcs
GRNCREPE120U	green	120 x 120	126 pcs
D WHCREPE050N	white	50 x 50	750 pcs
WHCREPE060N	white	60 x 60	750 pcs
WHCREPE075NU	white	75 x 75	348 pcs
WHCREPE090U	white	90 x 90	252 pcs
WHCREPE100U	white	100 x 100	252 pcs
WHCREPE120U	white	120 x 120	126 pcs
WHCREPE1315U	white	130 x 150	126 pcs



A



A TRAY LINERS

Our Tray Liners have high absorbency and rapid dispersion to help remove condensation which encourages dry packs and an effective sterilization process. To be used above or beneath trays. 80 gsm thickness. Compatible with steam sterilization. Single use.

REF	Format (mm)	Qty
TLC004	200 x 305	2000 pcs
TLC005	250 x 250	2400 pcs
TLC006	250 x 500	1200 pcs
TLC007	305 x 410	1000 pcs
TLC008	457 x 550	800 pcs
TLC009	305 x 305	2000 pcs
TLC010	305 x 457	1000 pcs
TLC011	305 x 610	1000 pcs
TLC012	457 x 457	1000 pcs
TLC013	457 x 610	600 pcs
TLC014	250 x 490	1200 pcs
TLC018	305 x 540	1000 pcs
TLC020	350 x 457	1000 pcs
TLC022	360 x 540	1000 pcs
TLC023	356 x 457	1000 pcs
TLC027	350 x 500	1000 pcs
TLC029	250 x 400	1200 pcs
TLC030	350 x 550	1000 pcs
TLC031	270 x 350	1000 pcs
TLC034	330 x 450	1000 pcs
TLC036	305 x 500	1000 pcs
TLC045	400 x 500	1000 pcs
TLC051	250 x 305	2000 pcs
TLC052	305 x 520	1000 pcs
TLC055	250 x 550	1200 pcs
TLC071	240 x 457	1000 pcs
TLC072	178 x 279	2400 pcs
TLC081	305 x 355	1000 pcs
TLC085	250 x 480	1000 pcs
TLC086	180 x 280	1000 pcs
TLC095	290 x 560	1000 pcs
TLC096	250 x 530	1200 pcs
TLC098	275 x 290	2000 pcs
TLC100	305 x 350	1000 pcs

A



A



A TRAY LINERS CONTINUED

Our Tray Liners have high absorbency and rapid dispersion to help remove condensation which encourages dry packs and an effective sterilization process. To be used above or beneath trays. 80 gsm thickness. Compatible with steam sterilization. Single use.

REF	Format (mm)	Qty
TLC102	160 x 260	2400 pcs
TLC108	350 x 610	1000 pcs
TLC109	130 x 230	2000 pcs
TLC110	150 x 250	2000 pcs
TLC114	250 x 415	1000 pcs
TLC115	315 x 525	1000 pcs
TLC116	457 x 500	1000 pcs
TLC122	457 x 630	400 pcs
TLC220	220 x 305	2000 pcs
TLC340	250 x 340	2000 pcs
TLC3535	350 x 350	1000 pcs
TLC3548	350 x 480	1000 pcs

A




Key Surgical tapes are designed for sealing sterilization packaging, such as fleece, paper, and crepe. They are offered with and without an indicator and are used to visually check whether the sterilization process has been carried out.

INDICATOR TAPES

For single use.

REF	Format (mm)	Adhesive	Color	Indicator	Color change	Qty
A 01506	25 x 50	strong	blue	steam	green to dark	1 roll
B 01500	19 x 50	strong	blue	steam	green to dark	1 roll
C 01505	24 x 50	standard	white	steam	green to dark	1 roll
D 01130	18 x 50	standard	white	steam (unleaded)	green to black	1 roll
E 01131	19 x 50	standard	beige	EO	red to green	1 roll
F 01151	19 x 50	standard	white	steam	green to brown	1 roll
G LF25TAPE 3390PLAS	25 x 25	standard	white	H ₂ O ₂	purple/blue to pink/ red	1 roll






ADHESIVE TAPES

For single use.

REF	Format (mm)	Adhesive	Color	Qty
A 01132	19 x 50	standard	white	1 roll
B 01507	25 x 50	standard	white	1 roll
C 01501	19 x 50	strong	blue	1 roll
D 01508	25 x 50	strong	blue	1 roll
E 01502	50 x 50	crepe	brown	1 roll





LABELED ADHESIVE TAPES


For single use, not heat resistant.

REF	Format (mm)	Caption	Font	Color	Qty
A 01527	50 x 66	INCOMPLETE	red	white	6 rolls
B 01798	50 x 50	USED MEDICAL INSTRUMENT	red	white	6 rolls
C 01799	50 x 50	DISINFECTED AND STERILE MEDICAL PRODUCTS	green	white	6 rolls

ACCESSORIES FOR ADHESIVE TAPES

REF	Description	Qty
D 09000	table dispenser for adhesive labels 120mm wide	1 pc
E 09001	table dispenser for 1 roll of indicator tape	1 pc
F 09002	table dispenser for 2 rolls of indicator tape	1 pc
09005	replacement knife for table dispenser REF 09001	1 pc
09007	replacement knife for table dispenser REF 09035	1 pc
09008	replacement knife for table dispenser REF 09002	1 pc
G 09035	table dispenser for indicator tape up to 25mm strip length can be preset (from approx. 1–10cm)	1 pc
H 09063	table dispenser for 4 rolls of indicator tape and tape up to 25mm, heavy version	1 pc
I 06013	table dispenser for adhesive labels for inside and outside rolls 250mm wide	1 pc

* Images in German only for reference



A



Wraps & Liners– Transport Wrap

K

Wraps & Liners

A TRANSPORT COVER BAGS

Our Transport Cover Bags are made of the same material as the Dry Transport Cover Wraps and can alternatively be sealed with indicator tape. Available in different colors (purple, red, yellow, and orange) and sizes.

REF	Format (mm)	Qty
TCOVER001	400 x 700 x 90	200 pcs
TCOVER002	500 x 700 x 90	200 pcs

B TRANSPORT WRAP DRY 60GSM

Protection against damaged packs, with strength and tear resistance with a patented absorbent layer. Particularly suitable for use when handling heavy orthopaedic sets. Available in purple, to distinguish from barrier wraps and also helps to support identification during transport and storage. Compatible with steam sterilization. Single use.

REF	Format (cm)	Qty
TDRY090	90 x 90	300 pcs
TDRY120	120 x 120	120 pcs

CHOICE BLUE AND TRANSPORT WRAP COMBINATIONS

REF	Format (cm)	Qty
CB/TDRY060	60 x 60	350 pcs
CB/TDRY090	90 x 90	300 pcs
CB/TDRY100	100 x 100	200 pcs
CB/TDRY120	120 x 120	120 pcs

B

