




# DISTAL PROTECTORS

The distal protection made of stretchable netted material adapts perfectly to the shape of nearly every endoscope. It is available in two sizes and is slipped onto the distal tip before transport to protect it from damage and costly repairs. The distal protection is sold non-sterile and individually packed for single use.

## ABOUT:

- Protection against damage to the endoscope distal tip
- Flexible fit for different endoscope sizes
- Easy to use



REF	DESCRIPTION	FORMAT	QTY
TP-100-03	Distal Protector to protect endoscope during transport, size: small, colour: red, non sterile	small	100 pcs.
TP-100-04	Distal Protector to protect endoscope during transport, size: large, colour: blue, non sterile	large	100 pcs.



TP-100-03 Distal protection in use with endoscope and 103614A-P transport tray

