



# ANTI-FOG SOLUTION

Don't let fog affect your sight during a surgical endoscopic procedure. Anti-Fog from Key Surgical is a perfect solution to help reduce condensation and fog on your laparoscopic and arthroscopic lens. The answer to a fog-free procedure has never been clearer.

## ABOUT:

- For surgical endoscopy procedures (e.g. laparoscopic and arthroscopic)
- Foam pad (radiopaque) has adhesive backing to fix on the sterile drape
- Sterile, single use and ready-to-use
- Non-toxic
- CE-marked - Class IIa
- Alcohol based solution (Isopropanol)



## FURTHER PRODUCT DETAILS

**REF AF-200**

6ml Anti-Fog solution, incl. 50 x 50mm foam pad

20 pcs/pack



a STERIS company

FOR FURTHER INFORMATION PLEASE CALL 49 4072374760 OR VISIT [KEYSURGICAL.DE](http://KEYSURGICAL.DE)

