



# POLYPECTOMY SNARES

Polypectomy Snares are used endoscopically in the removal of diminutive polyps, sessile polyps, pedunculated polyps and tissue from within the GI tract. Available in multiple styles (oval, mini-oval, or hexagonal) to suit your procedural needs.

## ABOUT:

- Self-contained metal loop that is opened over the polyp and then closed to capture and remove the polyp tissue
- Available in the most preferred styles: oval and hexagonal
- Braided wire combines the strength of a firm snare with the flexibility of a soft snare to help with accurately capturing tissue. This combination can also help reduce the need for additional snare inventory in the GI suite
- Sterile
- Single Use
- 10/box



PRODUCT CODE	STYLE	BIOPSY CHANNEL	WORKING LENGTH	SNARE DIAMETER	
EN10300	Stiff Multi-Wire, Hexagonal	2.8mm	2300mm	10mm	
EN10301	Stiff Multi-Wire, Hexagonal	2.8mm	2300mm	15mm	
EN10302	Stiff Multi-Wire, Hexagonal	2.8mm	2300mm	20mm	
EN10303	Stiff Multi-Wire, Hexagonal	2.8mm	2300mm	25mm	
EN10304	Stiff Multi-Wire, Hexagonal	2.8mm	2300mm	30mm	
EN10305	Stiff Multi-Wire, Hexagonal	2.8mm	2300mm	35mm	
EN10306	Stiff Multi-Wire, Oval	2.8mm	2300mm	10mm	
EN10307	Stiff Multi-Wire, Oval	2.8mm	2300mm	15mm	
EN10308	Stiff Multi-Wire, Oval	2.8mm	2300mm	20mm	
EN10309	Stiff Multi-Wire, Oval	2.8mm	2300mm	25mm	
EN10310	Stiff Multi-Wire, Oval	2.8mm	2300mm	30mm	
EN10311	Stiff Multi-Wire, Oval	2.8mm	2300mm	35mm	







FOR MORE INFORMATION CALL 01628 810626 OR VISIT [KEYSURGICAL.CO.UK](http://KEYSURGICAL.CO.UK)






# POLYPECTOMY SNARES


PRODUCT CODE	STYLE	BIOPSY CHANNEL	WORKING LENGTH	SNARE DIAMETER
EN10306R	Stiff Multi-Wire, Oval Rotatable	2.8mm	2300mm	10mm
EN10307R	Stiff Multi-Wire, Oval Rotatable	2.8mm	2300mm	15mm
EN10308R	Stiff Multi-Wire, Oval Rotatable	2.8mm	2300mm	20mm
EN10309R	Stiff Multi-Wire, Oval Rotatable	2.8mm	2300mm	25mm
EN10310R	Stiff Multi-Wire, Oval Rotatable	2.8mm	2300mm	30mm
EN10311R	Stiff Multi-Wire, Oval Rotatable	2.8mm	2300mm	35mm
EN10315	Multi-Wire, Crescent	2.8mm	2300mm	10mm
EN10316	Multi-Wire, Crescent	2.8mm	2300mm	15mm
EN10317	Multi-Wire, Crescent	2.8mm	2300mm	20mm
EN10318	Multi-Wire, Crescent	2.8mm	2300mm	25mm
EN10240	Monofilament Wire, Oval	2.8mm	2300mm	10mm
EN10241	Monofilament Wire, Oval	2.8mm	2300mm	15mm
EN10242	Monofilament Wire, Oval	2.8mm	2300mm	20mm
EN10243	Monofilament Wire, Oval	2.8mm	2300mm	25mm
EN10244	Monofilament Wire, Oval	2.8mm	2300mm	30mm
EN 10330R	Hybrid – Cold/Hot Rotatable	2.8mm	2300mm	10mm
EN 10331R	Hybrid – Cold/Hot Rotatable	2.8mm	2300mm	15mm



**YOU MIGHT ALSO LIKE**



Hemostasis Clips



Polyp Traps

