




ENDOSCOPE CLEANING BRUSHES

Endoscopes are complex and often difficult to reprocess so choosing the right brushes is critical to effective manual cleaning. Whether you need to clean a channel, a valve or a brush assortment, we got what you need, conveniently packaged and single use to reduce the risk of cross-contamination.

ABOUT:

- Available in a wide range of diameters and lengths, as a double or triple cleaning brush, for valves or the biopsy channel, single or as a convenient kit
- The ball headed tip of the brush protects the channel from perforations during cleaning
- Stable catheter with good grip
- Safe handling of wet and slippery instruments





ENDOSCOPE CLEANING BRUSHES

REF	DESCRIPTION	BRISTLE	WORKING LENGTH	LENGTH	QTY
EB10103	Double-Ended Brush	6.0mm	2.8mm	2500mm	50 pcs
EB10103.3	Triple-Cleaning Brush	6.0mm	2.8mm	2500mm	50 pcs
EB10136	Double-Ended Brush (Broncho)	3.0mm	1.4mm	1200mm	100 pcs
EB10102	Double-Ended Brush	3.0mm	2.0mm	2500mm	50 pcs
EB10104	Double-Ended Brush with Handle	5.0/10mm	2.8 - 4.2mm	2300mm	25 pcs
EB10100	Double-Ended Brush	5.0/10mm	2.8 - 4.2mm	2500mm	50 pcs

REF	ASSORTMENT	QTY
EB10107C	Brush Kit - Colonoscope (incl. Double-Ended Brush EB10103 & Channel Brush)	50 kit
EB10107G	Brush Kit - Gastroscope (incl. Double-Ended Brush EB10102 & Channel Brush)	50 kit
EB10103.3C	Brush Kit - Colonoscope (incl. Triple-Cleaning Brush EB10103.3 & Channel Brush)	50 kit

REF	DESCRIPTION	BRISTLE	WORKING LENGTH	LENGTH	QTY
EB10105	Pipe Cleaner, Double Ended	5.0mm	2,8 - 4.2mm	2300mm	50 pcs
EB10106	Cleaning Brush - Broncho	2.0mm	1.2mm	1200mm	50 pcs
*NEW EB10106.2L	Cleaning Brush - Broncho	2.0mm	1.2mm	2300mm	100 pcs
EB10108	Albarran Entrance Brush	2.3mm		13mm	100 pcs
26027	Disposable valve cleaning brush	5.2 / 11.0mm		156mm	100 pcs

