



STOPCON DISPOSABLE VALVES AND BIOPSY KITS

STOPCON Valves and Biopsy Kits for flexible endoscopes are provided sterile and ready-to-use. Available individually or conveniently packaged in a kit to provide you with exactly what you need, when you need it. The use of our STOPCON endoscope valves are single, but the benefits are many.

- No need for reprocessing/high-level disinfection. Single-use helps reduce cross-contamination.
- Oil-free valves (and no oil lubrication required). The use of silicone-based oil can create an environment for biofilm formation.

OLYMPUS®





Compatible with Olympus® Endoscopes 140/160/180/190/240140/160/180/190/240/260/290

	REF	DESCRIPTION
	EN10229.2	Air/Water and Suction Valves
	EN10229	Biopsy, Air/Water and Suction Valves
	EN10229.1	Biopsy, Air/Water, Suction Valves and Jet Channel Connector
*NEW 	EN10229.4	Biopsy, Air/Water, Suction Valves and Air/Water Cleaning Adapter
*NEW 	EN10229.5	Biopsy, Air/Water, Suction Valves, Jet Channel Connector and Air/Water Cleaning Adapter









OLYMPUS® AND FUJIJON®

	REF	DESCRIPTION
	EN10230	Air/Water Valve for Olympus®
	EN10232	Suction Valve for Olympus®
	EN10906	Jet Channel Connector for Olympus®
	EN10920	Air/Water Cleaning Adapter for Olympus®
 *NEW	EN 10933	Biopsy Port Valve for Olympus® Bronchoscopes
	EN10233	Biopsy Port Valve for Olympus® and Fujinon®
	EN10909	Biopsy Port Valve with Auxiliary Connector for Olympus® and Fujinon®
	EN10908	Jet Channel Connector for Fujinon®
	EN10229F	Biopsy, Air/Water and Suction Valves for Fujinon® 700

PENTAX®

	REF	DESCRIPTION
	EN 10234	Biopsy Port Valve
	EN10907	Jet Channel Connector
	EN10229P	Biopsy, Air/Water and Suction Valves
	EN10229P:1	Biopsy, Air/Water, Suction Valves and Jet Channel Connector

