




RELIANCE HEMOSTASIS CLIPS

Reliance Hemostasis Clips are used to prevent bleeding in the GI tract, which most often occurs after a polyp has been removed during a GI procedure. The clip is a small metal device that can be rotated and/or repositioned for accurate placement to join the surrounding tissue together.

ABOUT:

- Designed to clip tissue together during a GI procedure
- Repositionable and rotatable for improved accuracy, workflow, and performance
- Available in 11mm and 16mm opening widths
- Radiopaque
- MR Conditional
- Sterile
- Single-use



PRODUCT NO.	DESCRIPTION	QUANTITY
EN10400P.16	Hemostasis Clip RELIANCE, 230 cm length, 16 mm (open width)	10/box
EN10400P	Hemostasis Clip RELIANCE, 230 cm length, 11 mm (open width)	10/box

YOU MIGHT ALSO LIKE:

