



INSTRUMENT PROTECTORS

The Key Surgical instrument protectors made of paper and film not only protect the user from injuries, but also prevent the soft packaging from being damaged by the instrument. The paper side of the instrument protector ensures vapor penetration, while the foil side enables easy identification of the instrument. With a simple insertion procedure, the instrument can be stored safely and clearly in the strainer basket.



BENEFITS:

- Prevents damage to the sterile barrier
- Protects the user from injuries caused by the instrument
- Allows quick and safe identification of the instruments in the tray

	PROPERTIES	DESCRIPTION	FORMAT	UNITS
01455	Small instrument protector	Design to protect the instruments in the Steam sterilizer.	120 x 32 mm	100 pcs.
01457	Medium instrument protector	Design to protect the instruments in the Steam sterilizer.	160 x 40 mm	100 pcs.

* Usable length x max. Instrument width





KEYSURGICAL®